450 Series Waterproof Handheld Meter

Ion Selective Electrode (ISE) Operation Instructions

EUTECH
INSTRUMENTS
 Technology Made Easy ...


OAKTON[®]

Model:

Getting Started/Connections

After installing (2) AA batteries and/or connecting the optional 110/220 VAC power supply, connect the desired sensors to the corresponding ports.

12 mm and 16 mm probes can utilize the **Grip-Clip™** to attach one or more sensors to a beaker and to the instrument as needed. The stand can be extended as shown above or used for wall-mounting.

Keypad Functions

	Press once to power ON in the mode that was previously used. Press again to turn backlight on for one minute or off. Hold for 3 seconds to power OFF.
	Toggle between measurement and calibration modes. In SETUP mode, BACK serves to return to the previous menu option or setting.
	Confirm calibration values in CAL mode. Confirm selections in SETUP mode. Freeze or release the measured reading.
	Customize instrument settings and preferences. (See also Setup Programs)
	Toggle between available measurement types.
	Save measurement into memory. Increase value or scroll up in SETUP or manual calibration.
	Recall saved values from memory. Decrease value or scroll down in SETUP or manual calibration.
	Send output data to printer or computer.

Setup Programs

To access the settings below, press **SETUP**. Up/down arrows will display the available options. Press **ENTER** to accept the desired setting, or **BACK** to return to the previous option and/or exit.

Configuration Options

- Ready indicator **ON / OFF** / or Automatic **HOLD** when stable
- Choose **°Celsius** or **°Fahrenheit**

ISE Standard Options

- Choose standard group:
 - GRP.1** for 0.10, 1.00, 10.0, 100.0 & 1000 (5 points)
 - GRP.2** for 0.20, 2.00, 20.0, 200.0 & 2000 (5 points)
 - GRP.3** for 0.50, 5.00, 50.0 & 500.0 (4 points)
- Choose number of calibration points

Select Calibration Due Reminder

- Set number of days from **0-60** for desired parameter

View Calibration Data (for the parameter being measured)

- Press **ENTER** to view each point that is calibrated

View Electrode Data (for the parameter being measured)

- Press **ENTER** to view electrode data

System Settings

- Data Logging:
 - MANUAL** upon key press only
 - TIMED** interval. Choose (**SEC / MIN / HOUR**) interval.
- Automatic shut off after 10 minutes. Choose **ON** or **OFF**.
- Clock Settings:
 - Date: Choose **USA** (MM/DD/YYYY) or **Euro** (DD/MM/YYYY).
 - Time: Choose (**24HR** or **12HR**). If 12HR, choose **AM** or **PM**.
- Set Printer Type:
 - CSV** (Comma Separated Values) – best format for computer.
 - Printer** (Text) – best format for printer.
- Choose Manual (**MAN**) upon key press or **TIMED** interval. If timed, choose (**SEC / MIN / HOUR**).

Reset

- NO**. Exits from reset menu options without action.
- FACTORY RESET**. Returns all settings except date/time and ATC calibration to factory default values after **ENTER** is pressed then restarts meter.
- DATA RESET**. Erases data stored in memory while retaining other settings after **ENTER** is pressed.
- CALIBRATION RESET**. Erases non-ATC calibration data while retaining other settings after **ENTER** is pressed.

ISE Calibration

The pH 450 can measure ion concentration such as ammonia or fluoride when using an ion selective electrode (ISE) for the specific ion of interest. See **Setup** to select one of (3) ISE standard groups based on your intended standards and range. Use a minimum of 2 values for calibration and prepare the corresponding ion calibration solutions. For best results always begin with your lowest standard value, followed by the next lowest standard.

Refer to your Ion Selective Electrode instruction manual for details on conditioning, storage, maintenance, calibration standard preparation, Ionic Strength Adjustment, troubleshooting, etc. Each ISE is unique and requires care and operation that is specific to the electrode and ion of interest.

- Connect the ISE and press **MODE** as needed to view ppm. The primary display will show “- - -” without a stored ion calibration.
- Rinse the ISE with clean water and dip into your lowest calibration standard and stir. Press **CAL**.
- Use **▼/▲** keys to match the primary display to your desired ion calibration standard (0.10, 1.0, 10, etc). The secondary display is the corresponding mV.
- Press **ENTER** when the reading is stable to accept the calibration. **ERR** indicates that the calibration for the current point was not successful. This occurs when the slope (mV change between standards) is <15 mV/decade or >90 mV/decade. The primary display will show the next highest calibration standard. Use **▼/▲** keys to select a different ion calibration standard if needed.
- Rinse the ISE then dip into the corresponding calibration standard and stir. Repeat Steps 3 & 4 as needed. If the number of calibration points has been reached, the display will return to measurement mode. **Note:** The mV/decade slope value will be displayed in during measurement mode on the lower display upon a successful calibration. The value corresponds to the slope of the active reading.

mV Offset Adjustment

- While in mV measurement mode, dip the ORP and ATC sensors into a solution with a known mV value (i.e. Zobel, Light's, quinhydrone, or iodide/triiodide) and stir.

- When the “**READY**” indicator appears, use up/down arrows to adjust the primary reading to match the mV value at the measured temperature, then press **ENTER**. The meter allows an adjustable maximum value of ± 200 mV from the factory default mV value. When an offset has been stored successfully, RmV replaces mV.

Temperature Calibration/Manual ATC

- Press **CAL** from any measurement, then press **MODE**.
- Skip to step 3 for manual ATC, otherwise, dip the temperature sensor into a solution with a known accurate temperature. The upper display shows the active temperature while the lower display shows the factory default temperature without adjustment.
- Use up/down arrows to adjust the upper display. Press **ENTER** to accept the calibration temperature. The maximum adjustable value is ± 10 °C (or ± 18 °F) from factory default.

Error Messages

- ERR** “ERR” will appear when an error condition exists or the incorrect key is pressed. Common examples include:
- Pressing **ENTER** during calibration before the “**READY**” indicator appears. Wait for the “**READY**” indicator before pressing **ENTER**.
 - The calibration point was not successful. The allowable mV difference is 15 to 90 mV/decade.
 - UR** (Under Range) • **OR** (Over Range)

Intended Use, Maintenance & Precautions

These handheld meters use sensors to detect various parameters for water-based measurements. For routine maintenance disconnect the power cord or battery, then dust or wipe the display using a damp cloth. If necessary, warm water or a mild water based detergent can be used. Immediately remove any spilled substance from contact with the meter using the proper cleaning procedure for the type of spill.

- Do not use this equipment in potentially explosive atmospheres.
- Refer to the electrode instructions for use, storage and cleaning.
- Ensure that no liquid enters the instrument.
- Do not use any aggressive cleaning chemicals (solvents or similar agents).
- There are no user serviceable parts inside. Attempts to service internal parts may void the warranty.
- WARNING:** No modification of this equipment is allowed.

Instrument Operating Conditions	
Operating Ambient Temp.	5 to 45 °C
Operating Relative Humidity	5 to 85 %, non-condensing
Storage Temp.	-20 to +60 °C
Storage Relative Humidity	5 to 85 %, non-condensing
Pollution	Degree 2
Overvoltage	Category II
Weight	500 g
Size (L x W x H)	21.15 x 9.87 x 5.85 cm
Regulatory & Safety	CE, TUV 3-1, FCC Class A
Power Rating	DC Input: 9 VDC 1 A
Battery Requirement	2 x AA (LR6) 1.5 V batteries (replace batteries when battery sign blinks)
Vibration	Shipping/handling per ISTA #1A
Shock	Drop test in packaging per ISTA #1A
Enclosure (Designed To Meet)	IP67 (using rubber covers)
Universal Power Adapter Operating Conditions	
Operating Ambient Temp.	0 to 50 °C
Operating Relative Humidity	0 to 90 %, non-condensing
Storage Temp.	-20 to +75 °C
Storage Relative Humidity	0 to 90 %, non-condensing
Pollution	Degree 2
Overvoltage	Category II
Power Rating	I/P: 100 - 240 V, 50/60 Hz, 0.3A O/P: 9 VDC 1 A

Eutech Instruments Pte Ltd
 Blk 55, Ayer Rajah Crescent,
 #04-16/24, Singapore 139949
 Tel: (65) 6778-6876
 Fax: (65) 6778-0086
 eutech@thermofisher.com
 www.eutechinst.com

Oakton Instruments
 625 East Bunker Court,
 Vernon Hills, IL, 60061, USA
 Tel: 1-888-462-5866
 Fax: 1-847-247-2984
 info@4oakton.com
 www.4oakton.com

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EUTECH
INSTRUMENTS
Technology Made Easy ...

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