

SPECIFICATIONS

Range:	0 to 60 minutes
Resolution:	1 minute
Alarm:	Audible
Size:	3 $\frac{3}{4}$ H x 3 $\frac{3}{4}$ W x 2D inches
Weight:	3 $\frac{3}{4}$ ounces

OPERATION

Turn the knob clockwise to 60 and then turn it back to the desired time.

The bell will ring at the end of the selected time.

CONTROL COMPANY

12554 Galveston Rd. Suite B230

Webster, Texas 77598 USA

Ph. 281 482-1714 • Fax 281 482-9448

E-mail sales@control3.com

www.traceable.com

Control Company is ISO 9001:2008 Quality-
Certified by DNV and ISO/IEC 17025:2005
accredited as a Calibration Laboratory by A2LA.

© Control Company Cat. No. 1062

©2016 © Control Company. 92-1062-10 Rev. 3 052416

60-MINUTE INTERVAL TIMER INSTRUCTIONS