

Refractometers for all Breweries

For Brix scale **PAL-1**, For Plato scale **PAL-Plato**



PAL-1

Reasons to choose refractometer

Much Easier than Hydrometer

- Only a small amount of sample required!
- Measurement time, only 3 seconds!
- Rinse under running water
- Attractive design !
- Digital display - eliminates subjectivity



PAL-Plato



Specifications

	PAL-1 Cat.No.3810	PAL-Plato Cat.No.4590
Measurement range	: Brix 0.0 to 53.0%	Plato 0.0 to 30.0 °P
Measurement accuracy	: Brix $\pm 0.2\%$	Plato $\pm 0.2^\circ\text{P}$
Repeatability	: Brix $\pm 0.1\%$	Plato $\pm 0.1^\circ\text{P}$
Measurement temperature	: 10 to 75 °C (Automatic Temperature Compensation)	
Ambient temperature	: 10 to 40 °C	
International Protection Class	: IP65 Water resistant	
Dimensions and weight	: 55(W) × 31(D) × 109(H)mm, 100g (Main Unit only)	



Quick & Easy Cleanup!

- These models are basically designed to measure sweet wort.
For the measurement during and/or after fermentation, please understand the difference from densimeter/hydrometer before use.

ATAGO CO., LTD.

<http://www.atago.net/> overseas@atago.net

Headquarters: The Front Tower Shiba Koen, 23rd Floor
2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan
TEL : 81-3-3431-1943 FAX : 81-3-3431-1945



HACCP GMP GLP

ATAGO products comply with HACCP, GMP, and GLP system standards.

ATAGO U.S.A., Inc.

TEL : 1-425-637-2107

customerservice@atago-usa.com

ATAGO INDIA Instruments Pvt. Ltd.

TEL : 91-22-28544915, 40713232

customerservice@atago-india.com

ATAGO THAILAND Co., Ltd.

TEL : 66-21948727-9

customerservice@atago-thailand.com

ATAGO BRASIL Ltda.

TEL : 55 16 3913-8400

customerservice@atago-brasil.com

ATAGO ITALIA s.r.l.

TEL : 39 2 36557267

customerservice@atago-italia.com

ATAGO CHINA Guangzhou Co., Ltd.

TEL : 86 20-38106065

info@atago-china.com

ATAGO RUSSIA Ltd.

TEL : 7-812-777-96-96

info@atago-russia.com

* Specifications and appearance are subject to change without notice.

ENV.04 1503K Printed in Japan

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084