



novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682
Tel: (281) 359-8538

Fax: (866) 433-6684
Fax: (281) 359-0084

Mortar & Pestles

VEE GEE®, made by , Porcelain Mortar & Pestles 53121-Series, are manufactured with the highest standards. Mortars are glazed outside and on inside lip. No glaze on bottom or inside grinding surface. Pestle handle glazed to grinding surface. Used for mixing, crushing and grinding.



Specifications

Cat. No.	Mortar/ Pestle	Capacity (mL)	Mortar O.D. (mm)	Pestle Length (mm)	Pack Qty.
53121-1110	Mortar	70	63	NA	1 (Each)
53121-3100	Pestle	NA	NA	115	1 (Each)
53121-1220	Mortar	160	90	NA	1 (Each)
53121-3200	Pestle	NA	NA	135	1 (Each)
53121-1320	Mortar	400	125	NA	1 (Each)
53121-3400	Pestle	NA	NA	175	1 (Each)
53121-1500	Mortar	650	150	NA	1 (Each)
53121-3500	Pestle	NA	NA	210	1 (Each)



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • email: sales@veege.com

No.VGBRC080713/53121
Printed in the USA