

# first+care™ latex examination gloves



medical



emergency



veterinary



pharmacy



dental

Snapshots of some potential firstcare™ users\* at work:

- A pathologist examining a new tissue sample.
- A doctor performing an exam of a lump in a patient's arm.
- A dental lab technician making a new crown.
- A nursing home aid washing an infirm resident's head.

In healthcare, emergency work, scientific research, or other clinical settings, you need barrier protection that delivers the highest standards of reliability. Consistently. You need firstcare™.

Firstcare™ delivers the most reliable latex glove protection for people who absolutely need to avoid bodily fluids, pathogens, or contamination. Its outstanding tensile strength and resilience means a firstcare™ latex glove keeps its integrity even during the most physically demanding tasks. Low modulus and high elasticity allows for a true fit, giving hands the widest range of motion for intricate procedures while also minimizing hand fatigue. And the extremely low 1.5 AQL means you can be confident about its consistent protection, glove after glove. Whenever high clinical standards for safety are paramount, you'll be glad you have firstcare™ at hand, ready to protect.

\* Please consult your company policy to determine the appropriate glove for a given application.

**novatech**  
INTERNATIONAL

sales@novatech-usa.com

www.novatech-usa.com

Tel: (866) 433-6682

Tel: (281) 359-8538

Fax: (866) 433-6684

Fax: (281) 359-0084



**TDO:TESTED** Like all Spring-Fill Industries gloves, firstcare™ is the product of Test-Driven Optimization<sup>SM</sup> (TDO), our structured production and testing process to ensure our formulation of NRL provides you with maximum protection without compromise.

# firstcare™

## latex examination gloves



better grip

Unsurpassed Hand Dexterity, with Better Grip, and Form-fitted Comfort

The latex glove solution designed for people who want to maintain superior grip control plus barrier protection against bodily fluids, pathogens, or contamination.

Like all Spring-Fill Industries gloves, firstcare™ Latex is the product of Test-Driven Optimization<sup>SM</sup>, our structured production and testing process that ensures our formulation of NRL provides maximum protection without compromise. What's more, its extremely low 1.5 AQL means confidence about its consistent protection, glove after glove.

In healthcare, lab work, or any setting where keeping the finest control of instrumentation is top priority, you need barrier protection that delivers the highest standards of reliability and protection. You need firstcare™ Latex.



TDO: TESTED

### CHARACTERISTICS:

MATERIAL: PREMIUM NATURAL RUBBER LATEX  
TYPE: NON-STERILE  
USAGE: FOR SINGLE USE ONLY  
SHAPE: AMBIDEXTROUS

COLOR: NATURAL  
CUFF: INVERTED, BEADED CUFF  
STANDARDS: MEETS OR EXCEEDS FDA, EN, ASTM STANDARDS, CE

### TECHNICAL DATA:

GLOVE INTERIOR: POWDERED  
LENGTH (MM/IN):  
XS, S, M, L & XL: 240MM +/- 5MM  
THICKNESS:  
CUFF (MM/MILS): 0.091/3.58  
PALM (MM/MILS): 0.134/5.28  
FINGER (MM/MILS): 0.145/5.71  
TENSILE STRENGTH:  
BEFORE AGING (MPA): 18  
AFTER AGING (MPA): 14  
ELASTICITY/ELONGATION:  
BEFORE AGING (%): ≥650  
AFTER AGING (%): ≥500  
POWDER CONTENT (MG/G): <10  
PROTEIN LEVEL (μ G/G): <200

### TECHNICAL DATA:

GLOVE INTERIOR: POWDER-FREE  
LENGTH (MM/IN):  
XS, S, M, L & XL: 240MM +/- 5MM  
THICKNESS:  
CUFF (MM/MILS): 0.105/4.14  
PALM (MM/MILS): 0.134/5.28  
FINGER (MM/MILS): 0.145/5.71  
TENSILE STRENGTH:  
BEFORE AGING (MPA): 18  
AFTER AGING (MPA): 14  
ELASTICITY/ELONGATION:  
BEFORE AGING (%): ≥650  
AFTER AGING (%): ≥500  
POWDER CONTENT (MG/G): <2  
PROTEIN LEVEL (μ G/G): <50

### POWDERED

### ORDER INFORMATION

FIRSTCARE™  
5 MIL LATEX EXAMINATION GLOVES,  
NON-STERILE, AQL 1.5

### PACKAGING:

100 GLOVES/DISPENSER BOX (BY  
WEIGHT) 10 DISPENSER BOXES/CASE

### ORDER NUMBERS:

EXTRA SMALL: FCPT105A-XS  
SMALL: FCPT105A-S  
MEDIUM: FCPT105A-M  
LARGE: EXTRA FCPT105A-L  
LARGE: FCPT105A-XL

### POWDER-FREE

### ORDER INFORMATION

FIRSTCARE™  
5 MIL LATEX EXAMINATION GLOVES,  
NON-STERILE, AQL 1.5

### PACKAGING:

100 GLOVES/DISPENSER BOX (BY  
WEIGHT) 10 DISPENSER BOXES/CASE

### ORDER NUMBERS:

EXTRA SMALL: FCFT105A-XS  
SMALL: FCFT105A-S  
MEDIUM: FCFT105A-M  
LARGE: FCFT105A-L  
EXTRA LARGE: FCFT105A-XL

# The right solution for every job



## Jobguard™

**Jobguard™ Poly** Gloves

**Jobguard™ Vinyl** Multi-Purpose Gloves

**Jobguard™ Blue Vinyl** Multi-Purpose Gloves

**Jobguard™ Stretch Vinyl** Multi-Purpose Gloves

**Jobguard™ Nitrile** Multi-Purpose Gloves

**Jobguard™ XM Latex** Multi-Purpose Gloves

## first+care™

**firstcare™ Vinyl** Examination Gloves

**firstcare™ Stretch Vinyl** Examination Gloves

**firstcare™ Nitrile** Examination Gloves

**firstcare™ Latex** Examination Gloves



**providing comfort and peace of mind...  
on budget, on time**

**novatech**  
INTERNATIONAL

[sales@novatech-usa.com](mailto:sales@novatech-usa.com)

[www.novatech-usa.com](http://www.novatech-usa.com)

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084