

first+care™ vinyl examination gloves



Snapshots of some potential firstcare™ vinyl users* at work:

- A hospice worker helps an infirm patient out of bed.
- An ER orderly takes a patient's pulse and temperature.
- A nursing home aid bathes an immobile resident.
- A hospital housekeeping worker cleans a patients room.

In healthcare and related direct service work environments, where you may want to avoid latex but still need dependable manual barrier protection and a natural form factor for excellent dexterity, you need firstcare™ vinyl.

Delivering reliable protection for workers with hands-on tasks, firstcare™ vinyl provides an effective blend of strength and flexibility. Its formulation of polyvinyl chloride resin (pvc) means excellent barrier protection for most short term, lower risk procedures. Its ergonomic design ensures that tasks requiring precise gripping control can be carried out safely and successfully. And its thickness is precisely calibrated to yield superior protection with excellent touch sensitivity. For outstanding latex-free barrier protection needed for the majority of low risk, direct service tasks, firstcare™ vinyl provides protection without compromise.

* Please consult your company policy to determine the appropriate glove for a given application.



medical



emergency



veterinary



pharmacy



nursing home

novatech
INTERNATIONAL

sales@novatech-usa.com

www.novatech-usa.com

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084



TDO:TESTED Like every Spring-Fill Industries glove, firstcare™ Vinyl is manufactured using Test-Driven OptimizationSM, our system of exacting specifications, ensuring our formulation of PVC provides maximum protection without compromise.

first+care™

vinyl examination gloves



examination

Firstcare™ Vinyl is the glove solution designed for healthcare and related workers who want reliable protection and a greater degree of sensitivity without latex. Just like every Spring-Fill Industries glove, firstcare™ Vinyl is manufactured using Test-Driven OptimizationSM, our system of exacting specifications, ensuring our formulation of PVC provides the maximum protection without compromise.

And the smooth texture makes it perfect for tasks where touch sensitivity is needed. For latex-free, direct service environments that require a balance of barrier protection and tactile acuity, you need firstcare™ Vinyl.



CHARACTERISTICS:

MATERIAL: POLYVINYL CHLORIDE (100% LATEX FREE)
TYPE: NON-STERILE, AMBIDEXTROUS
USAGE: FOR SINGLE USE ONLY
CUFF: INVERTED, BEADED CUFF

COLOR: NATURAL/CLEAR*
STANDARDS: MEETS OR EXCEEDS FFDC AND FDA REQUIREMENTS, EN, ASTM STANDARDS, CE

TECHNICAL DATA:

GLOVE INTERIOR: POWDERED

LENGTH (MM/IN):
S, M, L & XL: 240mm +/- 5mm

THICKNESS:
CUFF (MM): 0.08 ± 0.02
PALM (MM): 0.09 ± 0.02
FINGER (MM): 0.10 ± 0.02

TENSILE STRENGTH:
BEFORE AGING (MPa): ≥11Mpa
AFTER AGING (MPa): ≥11Mpa

ELASTICITY/ELONGATION:
BEFORE AGING (%): MIN. 350
AFTER AGING (%): MIN. 350

POWDER CONTENT: MAX. 10mg/dm²
PROTEIN LEVEL: NONE

POWDERED

ORDER INFORMATION

FIRSTCARE™
VINYL EXAMINATION GLOVES
NON-STERILE

PACKAGING:

100 GLOVES/DISPENSER BOX (BY WEIGHT) 10 DISPENSER BOXES/MASTER CASE

ORDER NUMBERS:

SMALL: FCVPS105A-S
MEDIUM: FCVPS105A-M
LARGE: FCVPS105A-L
EXTRA LARGE: FCVPS105A-XL

TECHNICAL DATA:

GLOVE INTERIOR: POWDER-FREE

LENGTH (MM/IN):
S, M, L & XL: 240mm +/- 5mm

THICKNESS:
CUFF (MM): 0.08 ± 0.02
PALM (MM): 0.09 ± 0.02
FINGER (MM): 0.10 ± 0.02

TENSILE STRENGTH:
BEFORE AGING (MPa): ≥11Mpa
AFTER AGING (MPa): ≥11Mpa

ELASTICITY/ELONGATION:
BEFORE AGING (%): MIN. 350
AFTER AGING (%): MIN. 350

POWDER CONTENT: MAX. 2mg /glove
PROTEIN LEVEL: NONE

POWDER-FREE

ORDER INFORMATION

FIRSTCARE™
VINYL EXAMINATION GLOVES
NON-STERILE

PACKAGING:

100 GLOVES/DISPENSER BOX (BY WEIGHT) 10 DISPENSER BOXES/MASTER CASE

ORDER NUMBERS:

SMALL: FCVFS105A-S
MEDIUM: FCVFS105A-M
LARGE: FCVFS105A-L
EXTRA LARGE: FCVFS105A-XL

The right solution for every job



Jobguard™

Jobguard™ Poly Gloves

Jobguard™ Vinyl Multi-Purpose Gloves

Jobguard™ Blue Vinyl Multi-Purpose Gloves

Jobguard™ Stretch Vinyl Multi-Purpose Gloves

Jobguard™ Nitrile Multi-Purpose Gloves

Jobguard™ XM Latex Multi-Purpose Gloves

first+care™

firstcare™ Vinyl Examination Gloves

firstcare™ Stretch Vinyl Examination Gloves

firstcare™ Nitrile Examination Gloves

firstcare™ Latex Examination Gloves



**providing comfort and peace of mind...
on budget, on time**



novatech
INTERNATIONAL

sales@novatech-usa.com

www.novatech-usa.com

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084