

PORTABLE JAR MIXER

POSITIVE GEAR DRIVE

VARIABLE SPEED AGITATION

DIFFUSED LIGHTING

STURDY-COMPACT-SELF CONTAINED FOR FIELD OR LABORATORY USE



- A small gear head motor drives the stainless steel agitators, variable from 6 RPM to 250 RPM via a panel mounted potentiometer (optional digital readout of the paddle speed is available)
- Diffused lighting from a false bottom and controlled via an ON-OFF switch aids in evaluation of mix samples for clarity, floc particle size, and liquid-solid separation
- The front panel of the sturdy carrying case is removable for easy access to the agitators and controls
- The hinged lid stores pipettes and has a heavy duty handle for transporting the Jar Mixer
- The false bottom supports the 1000 ml mixing beakers and test solutions
- Includes a 120 VAC 60 Hz GFI protected power supply cord (optional 220 VAC 50/60 Hz available)
- The Portable Jar Mixer is available in two, four, and six station models

novatech
INTERNATIONAL

sales@novatech-usa.com

www.novatech-usa.com

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084