

Model	PMB 53	PMB 163	PMB 202
Capacity	50g	160g	200g
Readability	0.001g / 0.01%	0.001g / 0.01%	0.01g / 0.05%
Pan Size	3.9" ø / 100mm ø		
Weighing Units	g, %M, %S, ATRO %M, ATRO %S		
Stabilization Time (sec)	4		
Interface	RS-232, USB		
Memory Settings	Internal memory stores up to 99 test results and 49 programmable user settings External USB drive for limitless data storage of test results and settings		
Calibration	External calibration		
Display	Backlit LCD with 0.9"/ 24mm-high digits		
Power Supply	Power cord		
Operating Temperature	0° to 40°C / 32° to 104°F		
Housing	Die-cast aluminum		
Overall Dim.	9.8"x14.2"x7.3" / 250x360x185mm (wxdxh)		
Net Weight	13.2lb / 6kg		

Accessories

Item number	Description
104008036	Anti-vibration table
1070010636	Temperature calibration kit
301208188	In-use wet cover
3074010266	RS-232 cable (null-modem)
3074010267	USB cable
1120011156	ATP thermal printer
3126011263	ATP thermal printer paper
3126011281	ATP thermal printer paper (pack of 10)
307140001	Disposable sample pans (pack of 250)
700100046	Kensington-type lock and cable
600002028	AdamDU data collection program
700660290	Calibration certificate
700100027	F1 1g - 100g Calibration Weight Set Includes case to store weights
700100014	F1 50g Calibration Weight
700100012	F1 200g Calibration Weight

Key Accessories



ATP Thermal Printer

Easy-to-use features combined with quick set-up equals a truly versatile thermal printer. A single push-button feed key and three easy-to-read LED indicators make operation simple.



Temperature Calibration Kit

The temperature calibration kit comes in a cushioned carrying case for protection and easy transportation. The temperature measuring range is from -50°C to 300°C, with readings in Celsius or Fahrenheit.



Disposable Sample Pans

When weighing chemicals or other messy items, the use of disposable pans can help keep the permanent pan clean.

Supplied by

sales@novatech-usa.com
www.novatech-usa.com
 Tel: (866) 433-6682 Fax: (866) 433-6684
 Tel: (281) 359-8538 Fax: (281) 359-0084

www.adamequipment.com