

CALIBRATION REPORT

Digital Hand-held "Pocket" Refractometer

PAL-1

Serial No. **H414816**

	True value	Measurement value
Distilled water	Brix 0.0 %	0.0 %
Sugar solution	Brix 30.0 %	30.1 %
Sugar solution	Brix 50.0 %	50.0 %

Re-calibration is recommended every year to ensure the specified accuracy level.



CALIBRATION REPORT

Digital Hand-held "Pocket" Refractometer

PAL-2

Serial No. **H910025**

	True value	Measurement value
Sugar solution	Brix 50.0 %	50.0 %

Re-calibration is recommended every year to ensure the specified accuracy level.



CALIBRATION REPORT

Digital Hand-held "Pocket" Refractometer

PAL-3

Serial No. **H713795**

	True value	Measurement value
Distilled water	Brix 0.0 %	0.0 %
Sugar solution	Brix 30.0 %	30.0 %
Sugar solution	Brix 60.0 %	60.0 %

Re-calibration is recommended every year to ensure the specified accuracy level.

