

# BEAKERS

## Graduated Beakers

*"No-Drip" Spout for Easy Pouring*

Griffin low form polypropylene beakers are translucent, chemical resistant and will stack without sticking together. Raised numeric graduations show approximate volume ( $\pm 5\%$ ). All sizes are steam autoclavable at 121°C (250°F) and microwavable. **Per pack.**



CATALOG NO.	CAPACITY	PACKS/CASE	QNTY/PACK
F26210-0000	50ml	4	12
F26211-0000	100ml	4	12
F26212-0000	150ml	4	12
F26213-0000	250ml	6	6
F26214-0000	400ml	6	6
F26215-0000	600ml	6	4
F26216-0000	1000ml	4	3
F26217-0000	2000ml	6	1
F26219-0000	4000ml	6	1

## Beaker Specifications

CATALOG NO.	GRADUATION	HEIGHT	TOP O.D.	BOTTOM O.D.
F26210-0000	Every 10ml	48.3mm (1.9")	50.8mm (2")	40.6mm (1.6")
F26211-0000	25ml	66.0mm (2.6")	64.3mm (2.53")	50.8mm (2")
F26212-0000	25ml	76.2mm (3")	73.2mm (2.88")	58.4mm (2.3")
F26213-0000	50ml	88.9mm (3.5")	83.3mm (3.28")	66.8mm (2.63")
F26214-0000	100ml	111.1mm (4.37")	96.5mm (3.8")	74.7mm (2.94")
F26215-0000	100ml	127.0mm (5")	111.0mm (4.37")	83.8mm (3.3")
F26216-0000	200ml	152.0mm (6")	130.2mm (5.13")	101.6mm (4")
F26217-0000	500ml	184.2mm (7.25")	165.1mm (6.5")	133.4mm (5.25")
F26219-0000	1000ml	235.0mm (9.25")	212.6mm (8.37")	171.5mm (6.75")

## Weighted Beaker/Flask Holders


*Eliminate Drifting and Tipping in Water Baths*

Weighted holder grips beakers or flasks firmly in four flared stainless steel prongs. The Vikem® vinyl-coated lead ring secured beneath the holder keeps container immersed in water, provides a low center of gravity and allows water to circulate underneath the beaker. Can be used in boiling liquid up to 135°C(275°F) for 1 hour.



CATALOG NO.	FITS BEAKER	WEIGHT	DIAMETER
F18309-0100	100ml	165gm (0.36lb)	60mm (2.4")
F18309-0250	250ml	237gm (0.5lb)	80mm (3.3")

## Large Volume Beakers and Pitchers



5 and 10 Liter Capacity

Now Available with or without Handle!

These graduated beakers and pitchers feature flared rims and spouts for quick, easy pouring when large volumes are handled. Made of polypropylene, they offer good chemical resistance and are steam autoclavable at 121°C (250°F). 5 liter beakers and pitchers are graduated every 100ml. 10 liter beakers and pitchers are graduated every 200ml. **Per each.**



### BEAKERS (WITHOUT HANDLE)

CATALOG NO.	CAPACITY	BOTTOM O.D.	TOP O.D.	HEIGHT
F26219-0005	5 liter	16.5cm (6½")	21.0cm (8¼")	26.0cm (10¼")
F26219-0010	10 liter	22.2cm (8¾")	26.6cm (10½")	31.1cm (12½")

### PITCHERS (WITH HANDLE)


**NEW PRODUCT!**

CATALOG NO.	CAPACITY	BOTTOM O.D.	TOP O.D.	HEIGHT
F28994-0000	5 liter	16.50cm (6½")	21.0cm (8¼")	26.0cm (10¼")
F28995-0000	10 liter	22.2cm (8¾")	26.6cm (10½")	31.1cm (12½")

Scienceware® Graduated Pitchers Are Ideal for Pouring Large Volumes

## Clips for Beakers and Pitchers

Easy-to-use, polypropylene clips affix to the rim of a beaker or pitcher to hold tubing, dialysis bags or thermometers.



CATALOG NO.	DESCRIPTION	HOLDS	BAG/CASE PACKS
A. F18235-0000	Tubing Holder	Up to 9.5mm (¾") O.D.	6/6
B. F18237-0000	Dialysis Bag Holder	On rim up to ¼" thick	6/6
C. F18241-0000	Thermometer Holder	3.2mm (¼") to 11mm (½") dia.	3/12