

PITCHERS

Holdfast® Graduated Cylinders

Grip Handle for Steady Pouring

This one-piece polypropylene cylinder with handle provides a secure, solid grip for safer pouring. The "Total Contained, Total Delivered" double scale is easy-to-read and accurate to +/-5%. A hexagonal base prevents tipping and rolling. Handle opening is 8.9 x 2.9cm (3 x 1 1/8"). Chemically resistant and autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	GRADUATION	HEIGHT	I.D.	O.D.
F28461-1000	1000ml	Every 10ml	17.3"	2.3"	2.5
F28461-2000	2000ml	Every 20ml	19.3"	3.3"	3.5

See Chemical Resistance Chart on Pages 428-430

Tall Form Pitchers

Gracefully proportioned, these pitchers are ideal for general laboratory use and offer well designed spouts and large handles helpful for use with wet or dry ingredients. Both polypropylene and Polymethylpentene (TPX®) are chemical resistant and shatterproof. Metric gradations are molded on the outside wall. Steam Autoclavable at 121°C (250°F).



A. Graduated Pitchers - Polypropylene

Large selection of sizes. Translucent and heat resistant to 121°C (250°F).

Per each, cases listed below.



CATALOG NO.	CAPACITY	CASE
F28989-0000	250ml	24
F28990-0000	500ml	24
F28991-0000	1000ml	24
F28992-0000	2000ml	12
F28993-0000	3000ml	12
F28994-0000	5000ml (5L)	6
F28995-0000	10,000ml (10L)	1



B. Graduated Pitchers - TPX®

Crystal clear and heat resistant to 177°C (350°F). Per each, cases listed below.  

CATALOG NO.	CAPACITY	CASE
F28982-0000	500ml	24
F28983-0000	1000ml	24
F28984-0000	2000ml	12

Tall Form Graduated Pitchers - Additional Measurements

CAPACITY	GRADUATION	BOTTOM O.D.	TOP O.D.	OVERALL HEIGHT
250ml	Every 5ml	5.3cm (2")	6.2cm (2 ¹ / ₁₆ ")	11.9cm
500ml	10ml	7.4cm (2 ⁷ / ₁₆ ")	8.1cm (3 ¹ / ₁₆ ")	13.5cm
1000ml	10ml	9.4cm (3 ¹ / ₁₆ ")	10.7cm (4 ¹ / ₄ ")	16.4cm
2000ml	20ml	11.1cm (4 ³ / ₁₆ ")	12.7cm (5")	21.3cm
3000ml	50ml	12.2cm (4 ⁷ / ₁₆ ")	14.3cm (5 ⁵ / ₁₆ ")	26.4cm
5000ml	100ml	16.0cm (6 ¹ / ₁₆ ")	17.5cm (6 ⁷ / ₁₆ ")	27.1cm
10000ml	200ml	22.2cm (8 ¹ / ₁₆ ")	27.0cm (10 ⁵ / ₁₆ ")	32.0cm



Closed Spout Pouring Pitcher

Safe, Easy Pouring

Rigid, lightweight pouring pitcher has generous fill opening and extended spout for safe, easy pouring. Half-dome top prevents splattering and spillage. Spout measures 21mm (0.83") I.D. and 25mm (1") O.D. Made from polypropylene, it is steam autoclavable at 121°C (250°F). Graduated in liters and ml as well as U.S. and Imperial quarts and fractions. **Per each, 6 per case.**



CATALOG NO.	CAPACITY	HEIGHT	DIAMETER
F28742-0000	2 liter	30cm (12")	13cm (5 ¹ / ₄ ")

Use with Funnels, Large Industrial Tanks and Containers