



Infrared Micro-Sterilizer



The BactiZapper™ has an asbestos-free core element that uses infrared heating to produce 815°C (1500°F). It completely sterilizes microorganisms in 5-7 seconds. The BactiZapper uses well-proven technology to provide instant sterilization of platinum inoculating loops, needles, glass tube/pipette mouths and various metal & borosilicate glass instruments. The ceramic heating element is safely enclosed in a perforated stainless steel guard chamber, protecting users from accidental contact with the extreme temperatures of the core area. It is extremely compact, convenient and can be used in anaerobic and aerobic chambers.

- + **Complete infrared sterilization in 5-7 sec.**
- + **Safe, convenient and no spatter**
- + **Ideal for platinum loops, needles and tube mouths**

Technical Data:

<i>Max Temperature:</i>	815°C / 1500°F
<i>Sterilization Time:</i>	5 - 7 seconds
<i>Chamber Diameter:</i>	0.55 in. / 14mm
<i>Chamber Length:</i>	3.94 in. / 100mm
<i>External Dimensions:</i>	4 (W) x 6 (D) x 6.5 (H) in. 10.2 (W) x 15.2 (D) x 16.5 (H) cm
<i>Electrical Data:</i>	120V, 60 Hz / 230V, 50/60 Hz*
<i>Weight:</i>	2.9 lbs. / 1.3 kg

Ordering Information:

B1000*	BactiZapper Infrared Micro-Sterilizer
B1000-01	Loop holder attachment
B1000-03	Slide dryer attachment for 3 slides

* To order an item in 230 V, add (-E) to item number.

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084