

# **Stainless Steel Laboratory Jacks**

These premium lab jacks feature stainless steel top plates, slip proof rubber feet, and a screw operated mechanism that is easy to adjust

| STAINLESS STEEL LABORATORY JACKS |                  |               |              |             |
|----------------------------------|------------------|---------------|--------------|-------------|
| Code                             | Top Plate        |               | Dynamic Load | Static Load |
| Code                             | 10p Flate        | Max. Open Ht. | Dynamic Load | Static Load |
| SSJCT-10                         | 3 x 2.8 in       | 6 in          | 4.4 lb       | 33 lb       |
| SSJCT-20                         | 6.5 x 5 in       | 11 in         | 15.4 lb      | 66 lb       |
| SSJCT-30                         | 8 x 8 in         | 11 in         | 15.4 lb      | 66 lb       |
| SSJCT-40                         | 10 x 10 in       | 13 in         | 15.4 lb      | 66 lb       |
| SSJCT-50                         | 12 x 12 in       | 19 in         | 15.4 lb      | 66 lb       |
| SSJCT-60                         | 15.75 x 15.75 in | 19.7 in       | 15.4 lb      | 66 lb       |



## **Standard Laboratory Jacks**

These standard laboratory jacks feature steel construction and a screw opeated mechanism that is easy to adjust. They are light weight yet very stable and strong.

| STANDARD LABORATORY JACKS |            |               |              |             |
|---------------------------|------------|---------------|--------------|-------------|
| Code                      | Top Plate  | Max. Open Ht. | Dynamic Load | Static Load |
| CH0642A                   | 6 x 5 in   | 9.5 in        | 1 lb         | 55 lb       |
| CH0642B                   | 7.5 x 8 in | 11.8 in       | 1 lb         | 55 lb       |
| CH0642C                   | 6 x 6 in   | 9.8 in        | 2.2 lb       | 55 lb       |

## **Cast Iron Retort Base and Rod Sets**



These sets are made from premium qulity metalware and built to provide maximum stability and durability. The bases are cast iron with a high grade non toxic powdered coating. The stainless steel rods have precise threading that results in a snug and stable fit into the bases.

| HEAVY DUTY RETORT BASE AND ROD SETS |           |                      |
|-------------------------------------|-----------|----------------------|
| Code                                | Base      | Rod                  |
| CH0653ARD1                          | 4 X 6 in  | 19.5 in (0.4 in dia) |
| CH0653BRD2                          | 5 x 8 in  | 19.5 in (0.4 in dia) |
| CH0653D1RD3                         | 6 x 9 in  | 24 in (0.4 in dia)   |
| CH0653E1RD4                         | 6 x 11 in | 35.5 in (0.5 in dia) |



These sets are made from premium qulity metalware and built to provide maximum stability and durability. The bases are non ferrous cast iron with a high grade non toxic powdered coating and can be nested in one another for storage. The stainless steel rods have precise threading that results in a snug and stable fit into the bases.

| A SHAPE RETORT STAND AND ROD SETS |                    |         |
|-----------------------------------|--------------------|---------|
| Code                              | Base (Large Sides) | Rod     |
| STDRDAB01                         | 7 in               | 19.5 in |
| STDRDAB02                         | 8.7 in             | 23.6 in |
| STDRDAB03                         | 11 in              | 29.5 in |





# **Premium Laboratory Clamps**

- Solid zinc die-cast construction made from Zamac 3 Alloy
  + aluminum alloy
- Complete non-ferrous construction; rust resistant and durable in even the toughest scientific environments
- Wide range of options for various laboratory applications
- High grade non-toxic powdered coating
- For additional Clamps and specifications visit eiscoindustrial.com

|           | PREMIUM LABORATORY CLAMPS                                 |
|-----------|---|
| Code      | Description   |
| PRCMP01   | Premium Clamp, 3 Finger, Vinyl Coated                     |
| PRCMP01BH | Premium Clamp, 3 Finger, Vinyl Coated, W/ Bosshead        |
| PRCMP02   | Premium Clamp, 4 Prong, Cork Lined                        |
| PRCMP02BH | Premium Clamp, 4 Prong, Cork Lined, W/ Bosshead           |
| PRCMP03   | Premium Clamp, 4 Prong, Vinly Coated, On Rod              |
| PRCMP03BH | Premium Clamp, 4 Prong, Vinly Coated, W/ Bosshead         |
| PRCMP04   | Premium Clamp, 2 Prong, Vinyl Coated, 2x Adj              |
| PRCMP04BH | Premium Clamp, 2 Prong, Vinyl Coated, W/ Bosshead, 2x Adj |
| PRCMP05   | Premium Clamp, 2 Prong, Vinyl Coated                      |
| PRCMP05BH | Premium Clamp, 2 Prong, Vinyl Coated, W/ Bosshead         |
| PRCMP06   | Premium Clamp, 3 Prong, Vinyl Coated, 2x Adj              |
| PRCMP06BH | Premium Clamp, 3 Prong, Vinyl Coated, W/ Bosshead, 2x Adj |
| PRCMP07   | Premium Clamp, Long, 3 Finger, Vinyl Coated               |
| PRCMP07BH | Premium Clamp, Long, 3 Finger, Vinyl Coated, W/ Bosshead  |
| PRCMP08   | Premium Clamp, 2 Prong, Cork Lined                        |
| PRCMP08BH | Premium Clamp, 2 Prong, Cork Lined, W/ Bosshead           |
| PRCMP09   | Premium Clamp, 3 Finger, Cork Lined                       |
| PRCMP09BH | Premium Clamp, 3 Finger, Cork Lined, W/ Bosshead          |



### **Premium Bossheads**

#### PREMIUM BOSSHEAD - DIE CASTED

Eisco Industiral Series Right Angle Clamp Holder, made from high grade zinc alloy, finished in non-toxic powder coating with 2 stainless steel thumb screws.

- Threaded metal one piece clamping screws
- Able to hold rods up to 16 mm in diameter
- Holds rods orthogonal (90 degrees) to one another.
- All metal constructions to carry heavy loads and prevent torque damage

| Code   | Description                   |
|--------|-------------------------------|
| PRBH01 | Premium Bosshead - Die Casted |



#### PREMIUM BOSSHEAD - METAL

Eisco Industrial Series Mini All Metal Bosshead Clamp, for rods up to 15 mm.

- Threaded metal one piece clamping screws
- Able to hold rods up to 15 mm in diameter
- Great for when you need to connect rods and apparatus in tight places without sacrificing strength
- Holds rods orthogonal (90 degrees) to one another.
- All metal constructions to carry heavy loads and prevent torque damage

| Code   | Description              |
|--------|--------------------------|
| PRBH02 | Premium Bosshead - Metal |



#### PREMIUM BOSSHEAD - HEAVY DUTY

Eisco Industrial Series Heavy Duty Powder Coated Bosshead for rods up to 16 mm.

- Thicker than normal 7 mm threaded clamping screws
- Thicker casting to support more load
- Accepts up to 16 mm rods
- Made of high strength powder coated alloy
- All metal constructions to carry heavy loads and prevent torque damage

| Code   | Description                   |
|--------|-------------------------------|
| PRBH03 | Premium Bosshead - Heavy Duty |



#### PREMIUM BOSSHEAD - UNIVERSAL SWIVEL

Eisco Industrial Series Universal Swivel Bosshead for rods up to 16 mm, adjusts to any angle (dual 360° rotation of bossheads).

- Threaded metal one piece clamping screws
- Able to hold rods up to 16 mm in diameter
- Both bossheads can swivel 360 degrees
- Made of high strength powder coated alloy
- All metal constructions to carry heavy loads and prevent torque damage

| Code   | Description                         |
|--------|-------------------------------------|
| PRBH04 | Premium Bosshead - Universal Swivel |



#### PREMIUM BOSSHEAD - SQUARE HEAD

Eisco Industrial Series Classic Die Cast Alloy Bosshead, for rods up to 15 mm.

- Threaded metal one piece clamping screws
- Able to hold rods up to 15 mm in diameter
- Holds rods orthogonal (90 degrees) to one another.
- All metal constructions to carry heavy loads and prevent torque damage

| Code   | Description                    |
|--------|--------------------------------|
| PRBH05 | Premium Bosshead - Square Head |



# **Premium Metalware Sets**



| METALWARE SET - MTST1 |                          |
|-----------------------|--------------------------|
| Includes              | Item                     |
| CH0652B               | Base 20 X 12.5 cm, 800 g |
| CH0657B               | Rod 600mm                |
| CH0663A               | Bosshead                 |
| CH0682                | Clamp                    |
| CH0672D               | Retort ring 75mm O.D.    |



| METALWARE SET - MTST5 |                                |  |
|-----------------------|--------------------------------|--|
| Includes              | Item                           |  |
| CH0649B               | Base 20 X 12.5 cm, 1200 g      |  |
| CH0657B               | Rod 600mm                      |  |
| CH0672D               | Retort ring 75 mm O.D.         |  |
| CH0688                | Clamp - three prong with clamp |  |



| METALWARE SET - CH0652BR8 |                          |
|---------------------------|--------------------------|
| Includes                  | Item                     |
| CH0652B                   | Base 20 x 12.5 cm, 800 g |
| CH0659A                   | Rod 500 x 12 mm dia.     |
| CH0674B                   | Cast iron ring 75 mm     |
| CH0674C                   | Cast iron ring 100 mm    |
| CH0674D                   | Cast iron ring 125 mm    |

## **Premium Metalware Sets**

| CHEMISTRY LAB STARTER SET - MTST7 |                              |
|-----------------------------------|------------------------------|
| Includes                          | Item                         |
| CH0649C                           | Base 25 x 16 cm, 1800 g      |
| CH0676A                           | Three prong adjustable clamp |
| CH0659C                           | Stainless steel rod 750 mm   |
| CH0665A                           | Bosshead                     |
| CH0674B                           | Cast ring of 75 mm dia       |
| CH0993-NG                         | Burner bunsen - Natural Gas  |
| CH0258B                           | Double burette clamp         |
| CH0645-SPL                        | Tripod stand                 |
| CH0674C                           | Cast iron ring of 100 mm dia |
| CH0100E                           | Burner tubing - 60 cm        |
| CH0675B                           | Two prong adjustable clamp   |
| CH0525B                           | Wire Guage 150 mmmm O.D.     |



|                             | LAD FOOENTIALO KIT. LADKIT |  |
|-----------------------------|----------------------------|--|
| LAB ESSENTIALS KIT - LABKIT |                            |  |
| Includes                    | Item                       |  |
| CH0649B                     | Base 20 x 12.5 cm, 1200 g  |  |
| CH0672D                     | Retort ring 75 mm O.D.     |  |
| CH0657B                     | Rod 600 mm                 |  |
| CH0258B                     | Burette clamp double       |  |
| CH0729B                     | Crucible tong              |  |
| CH0713A                     | Test tube stand Z shape    |  |
| BKTG                        | Beaker tong                |  |
| CH0730                      | Test tube holder           |  |
| CH0644-SPL                  | Tripod stand               |  |
| CH0204G                     | Test tube brush            |  |
| CH0525A                     | Wire gauge                 |  |
| CH0632B                     | Spatula                    |  |
| CH0663A                     | Bosshead                   |  |
| CH0682                      | Clamp retort - 4 prong     |  |



#### ULTIMATE LAB STARTER KIT - ULTLABKIT

Includes the following items:

- 20 X 12.5cm Powder Coated Metal Base with Rubber Feet
- 600mm Stainless Steel Rod (Threaded for Base)
- Beaker Tongs
- Crucible Tongs
- Tripod Stand
- Bosshead
- Test Tube Clamp
- Sheet Metal Test Tube Stand

- Cork Lined Test Tube Holder
- Test Tube Brush
- Stainless Steel Spatula
- 2 Polypropylene 100ml graduated cylinders
- 2 Polypropylene 500ml beakers
- 50ml Acrylic Burette
- 250 ml borosilicate Erlenmeyer flask
- 14 X 12 X 11 cm cleaning/drying basket





Tel.: 585.563.7963 Fax.: 585.486.1735

E-mail: customer.support@eiscolabs.com