

- **VersaCap®** technology provides the flexibility of interchangeable connectors
- **EZgrip™** body features advanced ergonomic design
- **EZfil™** wide-mouth design for easy filling and cleaning
- **VersaBarb®** technology provides interchangeable spigot options

- **Made in the USA**



Caution

- Do not use carboys under pressure or full vacuum
- Do not use carboys at temperatures greater than the rated maximum
- Do not store oxidizing agents in carboys
- Do not flame carboy
- Do not use carboy spigot as a handle
- Do not use for long-term storage of strong organic compounds

Foxx Life Sciences' **EZgrip™** Carboy is an advanced fluid handling system developed with the end-user in mind. The ergonomically designed **EZgrip™** body, coupled with our **VersaCap®** technology provides the user with an easy-to-use and easy-to-clean carboy.

The **EZgrip™** HDPE Carboy is made from a translucent USP Class VI, FDA Grade resin, with excellent chemical resistance, making it ideal for most buffers and reagents as well as powders and solids.


The **EZgrip™** Carboy's body includes design features that allow the user to easily hold and manipulate during lifting, filling, and pouring, while **VersaCap®** technology provides interchangeable, leak-proof top connections. **EZgrip™** Carboy features include large top handles, molded hand grips in the base, easy-to-read metric graduation marks certified to $\pm 5\%$ accuracy, and the carboy material molded into the side for easy identification.

The **EZgrip™** HDPE Carboy offers an **EZfil™** 120mm neck for easier filling and cleaning. This is available on 20L, 40L, and 75L models.



- 1 Novel shape maximizes bench and shelf space
- 2 Bottom handles provide a secure grip while lifting and transporting
- 3 Large top handles for a full hand grip
- 4 Molded carboy identification marks
- 5 Large, clear graduation marks with $\pm 5\%$ accuracy
- 6 Autoclavable spigot and **VersaBarb®** technology
- 7 **VersaCap®** technology with versatile adapters and replaceable gasket
- 8 **EZgrip™** body design allows for maximum solution recovery
- 9 **EZfil™** wide-mouth shape for easy filling and cleaning

ISO 13485



Product Number	Description	Capacity (L)	Material	Spigot	Cap Size (mm)	Qty/Pack
152-0111-OEM	Carboy, HDPE, 2.5L, VersaCap® 83B	2.5	HDPE	No	83	1
152-1111-OEM	Carboy, HDPE, 5L, VersaCap® 83B	5	HDPE	No	83	1
152-1211-OEM	Carboy, HDPE, 5L, Spigot, VersaCap® 83B	5	HDPE	Yes	83	1
152-2111-OEM	Carboy, HDPE, 10L, VersaCap® 83B	10	HDPE	No	83	1
152-2211-OEM	Carboy, HDPE, 10L, Spigot, VersaCap® 83B	10	HDPE	Yes	83	1
152-3111-OEM	Carboy, HDPE, 20L, VersaCap® 83B	20	HDPE	No	83	1
152-3211-OEM	Carboy, HDPE, 20L, Spigot, VersaCap® 83B	20	HDPE	Yes	83	1
152-3121-OEM	Carboy, HDPE, 20L, VersaCap® 120mm	20	HDPE	No	120	1
152-3221-OEM	Carboy, HDPE, 20L, Spigot VersaCap® 120mm	20	HDPE	Yes	120	1
152-5121-OEM	Carboy, HDPE, 40L, VersaCap® 120mm	40	HDPE	No	120	1
152-5221-OEM	Carboy, HDPE, 40L, Spigot VersaCap® 120mm	40	HDPE	Yes	120	1
152-7121-OEM	Carboy, HDPE, 75L, No Spigot VersaCap® 120mm	75	HDPE	No	120	1



novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084