

Haier Biomedical Liquid Nitrogen Storage Solutions

Biobank Series LN₂ Solutions



- Large scale capacity stores up to 94,875×2ml Vials
- Vapour phase storage at -190°C prevent cross contamination
- Supplied with intelligent Cryosmart liquid nitrogen system monitoring and control solution
- LN₂ splash proof design for safer operation

Smart Series LN₂ Solutions



- 5 models from 2,400 to 6,000 cryovial capacity
- Cloud data storage for traceability
- Low consumption rate and high performance stability
- Durable aluminium construction
- Real time monitor of temperature and fluid level
- New lock design
- 5-year vacuum warranty

Medium Series LN₂ Solutions



- Heavy duty lockable enclosure
- Compatible with all major cryobox brands
- Liquid or vapour phase option available
- Durable aluminium construction
- Temperature monitoring available
- High thermal efficiency
- Ultra-low evaporation loss
- 5-year vacuum warranty

Dryshipper Series LN₂ Solutions



- Rugged and durable aluminium construction
- Vapour phase Storage
- Straw and cryovial storage options
- No LN₂ spillage
- Lockable lids
- Cryo absorbant
- 3-year Vacuum Warranty

Self-pressurized Series LN₂ Solutions



- Stainless steel construction
- Integral safety mechanisms
- Labelled valves for easy identification
- Safety valves
- Decant valves
- Pressure raising
- Liquid storage
- Rotary ring construction
- 5-year vacuum warranty

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International