

### Maximize Safety in your Lab

Disposable, highly absorbent paper with polyethylene backing is designed to protect laboratory surfaces against liquid spills and prevent leakage through to the bench top, necessary to maintain a safer work environment.



### FEATURES

- Meets OSHA Regulation 29CFR 1910.1030 for Occupational Exposure to Bloodborne Pathogens
- Tear resistant material, whether wet or dry
- Skid resistant, waterproof, chemical resistant polyethylene backing
- Smooth white surface for writing on
- Spilled or drained fluids are quickly trapped in absorbent paper
- Available in rolls or sheets

### APPLICATIONS

- Surface protector
- Contain spillage of hazardous substances to avoid contamination
- Water or solvent wick for humidity chambers and fume hoods
- Shelf and drawer liner
- Economical way to extend utility or maintenance carts



sales@novatech-usa.com  
www.novatech-usa.com  
Tel: (866) 433-6682 Fax: (866) 433-6684  
Tel: (281) 359-8538 Fax: (281) 359-0084

# Bench Protector Paper

## TECHNICAL SPECIFICATIONS:

- Material:** Pure cellulose with polyethylene backing  
**Basis weight:** 140 g/m<sup>2</sup> (86 lbs/3,000 ft<sup>2</sup>) – TAPPI T410  
**Thickness:** 0.22 mm (0.0086 in, 8.66 mil) – TAPPI T411  
**Surface:** (A) Smooth, bright white highly absorptive  
(B) Smooth polyethylene clear laminate

## ORDERING INFORMATION

| Tremont Cat#      | Description   |
|-------------------|---|
| <b>LEBPR4850M</b> | Bench protector paper roll 48 cm x 100 m, 1 roll per pack       |
| <b>LEBPS4860</b>  | Bench protector paper sheets 48 cm x 60 cm, 100 sheets per pack |

## COMMON COMPETITIVE EQUIVALENTS

| Dimensions           | Tremont Cat#      | Company  | Catalog#  |
|----------------------|-------------------|--|-----------|
| 48 cm x 100 m Roll   | <b>LEBPR4850M</b> | Kimberly-Clark Kimtech Science®<br>Bench Top Protector | 75450     |
|                      |                   | GE Healthcare/Whatman Benchkote™                       | 2300-731  |
|                      |                   | Ahlstrom Polyshield®                                   | 4002-4650 |
| 48 cm x 60 cm Sheets | <b>LEBPS4860</b>  | Kimberly-Clark Kimtech Science®<br>Bench Top Protector | 75460     |
|                      |                   | GE Healthcare/Whatman Benchkote™                       | 2300-917  |
|                      |                   | Ahlstrom Polyshield™                                   | 4002-4657 |

## PACKAGING

Bench Protector Paper is packaged in roll or sheet formats.



©2013 I.W. Tremont Co. Inc.  
ISO9001:2008 registration by:



Trademarks are owned by their respective owners. Use of competitive trademarks are for reference only – no ownership or rights are implied.

Cross reference information is a compilation of the listed manufacturers equivalency charts and sales data, not actual test data.

Although reasonable attempt has been made to ensure equivalency between competitive products - individual testing is suggested for non-standard method applications.

**I.W. Tremont**

**Technical and Specialty Papers**

ISO 9001:2008 Registered