

Petri Dish

I.W. Tremont model #: LEMB-2050 – Petri Dishes LEMB-2051 – Petri Dishes with cellulose pads

Description:

	LEMB-2050	LEMB-2051
1. Diameter	55mm	55mm
2. High	15mm	15mm
3. Material	Crystal polystyrene	Crystal polystyrene
4. Cellulose Pad	No	Yes
5. Sterility Mode	Gamma radiation	Gamma radiation
6. Quantity per box	720 units	100 units

Known applications:

Used to culture bacteria, fungi, and other microorganisms. Petri dish with absorbent pad eliminates hand-loading time and minimizes contamination.

Closest equivalents:

I.W. Tremont#	LEMB-2050		LEMB-2051	
	Part #	Same / (difference)	Part #	Same / (difference)
Company:				
Millipore	PD2004705	(47mm diameter, 600/pack)	PD20047S0	(47mm diameter, 150/pack)
Pall	7232	(50mm diameter, 500/pack)	7245	(50mm diameter)
Advantec	800500	(50mm diameter, 500/pack)	800101	(50mm diameter)
Fisher	09-720-502	(47mm diameter, 600/pack)	09-720-501	(47mm diameter, 150/pack)
VWR	25-384-092	(60mm diameter, 500/pack)		
Nunc	150326	(60mm diameter, 400/pack)		
Corning/Gosselin	BP53-02	(60mm diameter, 1620/pack)		