

## Boston Round Glass Bottles



- Available in clear or amber borosilicate glass

- Choice of closure and liner

- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

**CLASS 1: (Standard)** Containers are assembled with liner and closure without washing treatment.

**CLASS 2: (Precleaned)** Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

**CLASS 3: (Certified)** Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

### Assembled with White Polypropylene Closure, .015" PTFE Lined

#### Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
9-169	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread, White Closure, PTFE Lined	24 Pack
9-170	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-171	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-172	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9-173	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9-174	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9-175	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

#### Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
9A-169	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-170	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-171	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-172	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-173	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-174	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-175	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

## Boston Round Glass Bottles

### Precleaned

Assembled with White Polypropylene Closure, .015" PTFE Lined

#### Clear Glass Boston Round Bottles - Class 2

Cat. No.	Description	Qty
9-169-2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread, White Closure, PTFE Lined	24 Pack
9-170-2	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-171-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-172-2	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9-173-2	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9-174-2	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9-175-2	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

Assembled with White Polypropylene Closure, .015" PTFE Lined

#### Amber Glass Boston Round Bottles - Class 2

Cat. No.	Description	Qty
9A-169-2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-170-2	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-171-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-172-2	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-173-2	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-174-2	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-175-2	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

### Precleaned and Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

#### Clear Glass Boston Round Bottles - Class 3

Cat. No.	Description	Qty
9-170-3	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-171-3	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-172-3	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9-173-3	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9-174-3	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9-175-3	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

J.G. Finneran Associates, Inc. — Visit [www.JGFinneran.com](http://www.JGFinneran.com) • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • [sales@novatech-usa.com](mailto:sales@novatech-usa.com) • [www.novatech-usa.com](http://www.novatech-usa.com)

## Boston Round Glass Bottles

Precleaned and Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

### Amber Glass Boston Round Bottles - Class 3

Cat. No.	Description	Qty
9A-170-3	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-171-3	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-172-3	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-173-3	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-174-3	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-175-3	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

**J.G. Finneran Associates, Inc.** — Visit [www.JGFinneran.com](http://www.JGFinneran.com) • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • [sales@novatech-usa.com](mailto:sales@novatech-usa.com) • [www.novatech-usa.com](http://www.novatech-usa.com)

## Boston Round Glass Bottles with Assorted Closures

### Standard

#### Assembled with Black Phenolic Closure/PE Cone Lined

##### Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0149-1/2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-4	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-8	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0149-16	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0149-32	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack

##### Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0161-1/2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-4	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-8	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0161-16	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0161-32	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack

#### Assembled with White Polypropylene Closure, F217 Lined

##### Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0156-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0156-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0156-4	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0156-8	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0156-16	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0156-32	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

## Boston Round Glass Bottles with Assorted Closures

### Standard

#### Assembled with White Polypropylene Closure, F217 Lined

##### Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0168-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0168-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0168-4	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0168-8	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0168-16	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0168-32	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

#### Assembled with White Polypropylene Closure, Unlined

##### Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0154-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0154-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0154-4	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0154-8	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0154-16	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0154-32	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

##### Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0166-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0166-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0166-4	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0166-8	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0166-16	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0166-32	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

**J.G. Finneran Associates, Inc.** — Visit [www.JGFinneran.com](http://www.JGFinneran.com) • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • [sales@novatech-usa.com](mailto:sales@novatech-usa.com) • [www.novatech-usa.com](http://www.novatech-usa.com)

## Boston Round Glass Bottles with Assorted Closures

### Standard

#### Without Closures

##### Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0155-1/2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread	24 Pack
D0155B-1/2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread - Bulk Pack	720 Pieces
D0155-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread	24 Pack
D0155B-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread - Bulk Pack	432 Pieces
D0155-2	2 oz, 60mL Clear Boston Round Bottle, Cap, 20-400mm Thread	24 Pack
D0155B-2	2 oz, 60mL Clear Boston Round Bottle, Cap, 20-400mm Thread - Bulk Pack	288 Pieces
D0155-4	4 oz, 125mL Clear Boston Round Bottle, Cap, 22-400mm Thread	24 Pack
D0155B-4	4 oz, 125mL Clear Boston Round Bottle, Cap, 22-400mm Thread - Bulk Pack	160 Pieces
D0155-8	8 oz, 250mL Clear Boston Round Bottle, Cap, 24-400mm Thread	12 Pack
D0155-16	16 oz, 500mL Clear Boston Round Bottle, Cap, 28-400mm Thread	12 Pack
D0155-32	32 oz, 1000mL Clear Boston Round Bottle, Cap, 33-400mm Thread	12 Pack

##### Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0167-1/2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread	24 Pack
D0167B-1/2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread - Bulk Pack	720 Pieces
D0167-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread	24 Pack
D0167B-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread - Bulk Pack	432 Pieces
D0167-2	2 oz, 60mL Amber Boston Round Bottle, Cap, 20-400mm Thread	24 Pack
D0167B-2	2 oz, 60mL Amber Boston Round Bottle, Cap, 20-400mm Thread - Bulk Pack	288 Pieces
D0167-4	4 oz, 125mL Amber Boston Round Bottle, Cap, 22-400mm Thread	24 Pack
D0167B-4	4 oz, 125mL Amber Boston Round Bottle, Cap, 22-400mm Thread - Bulk Pack	160 Pieces
D0167-8	8 oz, 250mL Amber Boston Round Bottle, Cap, 24-400mm Thread	12 Pack
D0167-16	16 oz, 500mL Amber Boston Round Bottle, Cap, 28-400mm Thread	12 Pack
D0167-32	32 oz, 1000mL Amber Boston Round Bottle, Cap, 33-400mm Thread	12 Pack

