

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall



- Ideal for use in soil, sediment and sludge sampling

- Preassembled with .015" PTFE or F217 lined closed top white polypropylene closures

- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
9-180	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-181	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-182	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-183	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-184	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
9-190	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-191	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-192	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-193	64 oz, 2000mL Clear Tall Jar, 83-400mm Thread, White Polypropylene Closure, PTFE Lined	6 Pack
9-194	128 oz, 4000mL Clear Tall Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall

Precleaned

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 2

Cat. No.	Description	Qty
9-180-2	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-181-2	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-182-2	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-183-2	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-184-2	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 2

Cat. No.	Description	Qty
9-190-2	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-191-2	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-192-2	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-193-2	64 oz, 2000mL Clear Tall Jar, 83-400mm Thread, White Polypropylene Closure, PTFE Lined	6 Pack
9-194-2	128 oz, 4000mL Clear Tall Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Precleaned and Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 3

Cat. No.	Description	Qty
9-180-3	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-181-3	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-182-3	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-183-3	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-184-3	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 3

Cat. No.	Description	Qty
9-190-3	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-191-3	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-192-3	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-193-3	64 oz, 2000mL Clear Tall Jar, 83-400mm Thread, White Polypropylene Closure, PTFE Lined	6 Pack
9-194-3	128 oz, 4000mL Clear Tall Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

J.G. Finneran Associates, Inc. — Visit www.JGFinneran.com • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • sales@novatech-usa.com • www.novatech-usa.com

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall with Assorted Closures

Standard

Assembled with White Polypropylene Closure, F217 Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0086-2	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0086-4	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0086-8	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0086-16	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0086-32	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0098-1	1 oz, 30mL Clear Tall Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0098-2	2 oz, 60mL Clear Tall Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0098-4	4 oz, 250ml Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0098-8	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0098-16	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

Assembled with White Polypropylene Closure, Unlined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0084-2	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0084-4	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0084-8	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0084-16	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0084-32	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0096-1	1 oz, 30mL Clear Tall Jar, 33-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0096-2	2 oz, 60mL Clear Tall Jar, 38-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0096-4	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0096-8	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0096-16	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall with Assorted Closures

Standard

Without Closures

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0085-2	2 oz, 60mL Clear Jar, 53-400mm Thread	24 Pack
D0085-4	4 oz, 125mL Clear Jar, 58-400mm Thread	24 Pack
D0085-8	8 oz, 250mL Clear Jar, 70-400mm Thread	24 Pack
D0085-16	16 oz, 500mL Clear Jar, 89-400mm Thread	12 Pack
D0085-32	32 oz, 1000mL Clear Jar, 89-400mm Thread	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0097-1	1 oz, 30mL Clear Tall Jar, 33-400mm Thread	24 Pack
D0097-2	2 oz, 60mL Clear Tall Jar, 38-400mm Thread	24 Pack
D0097-4	4 oz, 125mL Clear Tall Jar, 48-400mm Thread	24 Pack
D0097-8	8 oz, 250mL Clear Tall Jar, 58-400mm Thread	12 Pack
D0097-16	16 oz, 500mL Clear Tall Jar, 70-400mm Thread	12 Pack

