

High Density Polyethylene Wide Mouth Jars



HDPE Wide Mouth Jars

- Preassembled with foamed polyethylene lined polypropylene closures
- All HDPE bottles and jars are also available with PTFE-lined closures or unlined closures.
- Available in standard, and precleaned in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

Standard

Assembled with White Polypropylene Closure, F217 Lined

HDPE Wide Mouth Jars - Class 1		
Cat. No.	Description	Qty
9-200	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-200B	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	960 Pieces
9-201	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-201B	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	540 Pieces
9-202	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-202B	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-203	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-203B	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-204	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-204B	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	72 Pieces
9-205	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined	6 Pack
9-205B	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	30 Pieces
9-206	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined	4 Pack
9-206B	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	12 Pieces

High Density Polyethylene Wide Mouth Jars

Precleaned

Assembled with White Polypropylene Closure, F217 Lined

HDPE Wide Mouth Jars - Class 2		
Cat. No.	Description	Qty
9-200-2	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-200B-2	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	960 Pieces
9-201-2	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-201B-2	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	540 Pieces
9-202-2	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-202B-2	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-203-2	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-203B-2	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-204-2	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-204B-2	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	72 Pieces
9-205-2	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined	6 Pack
9-205B-2	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	30 Pieces
9-206-2	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined	4 Pack
9-206B-2	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	12 Pieces

J.G. Finneran Associates, Inc. — Visit www.JGFinneran.com • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • sales@novatech-usa.com • www.novatech-usa.com

High Density Polyethylene Narrow Mouth Bottles



HDPE Narrow Mouth Bottles

- *Preassembled with foamed polyethylene lined polypropylene closures*
- *All HDPE bottles and jars are also available with PTFE - lined closures or unlined closures.*
- *Available in standard, and precleaned in accordance with recommended E.P.A. Protocol*

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Pre-cleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

Standard

Assembled with White Polypropylene Closure, F217 Lined

HDPE Narrow Mouth Bottles - Class 1

Cat. No.	Description	Qty
9-210	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-210B	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	470 Pieces
9-211	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-211B	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-212	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-212B	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-213	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-213B	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	84 Pieces

Precleaned

Assembled with White Polypropylene Closure, F217 Lined

HDPE Narrow Mouth Bottles - Class 2

Cat. No.	Description	Qty
9-210-2	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-210B-2	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	470 Pieces
9-211-2	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-211B-2	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-212-2	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-212B-2	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-213-2	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-213B-2	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	84 Pieces