

Laboratory Grade Plastic Bottles



- Available in wide mouth or narrow mouth

- Available in natural or amber

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

Standard

Wide Mouth Plastic Laboratory Grade Bottles with Closures

High Density Polyethylene (HDPE), Natural - Class 1

Cat. No.	Description	Qty
D0433-4	125mL Natural Wide Mouth Laboratory Grade Bottle, 38mm Closure	24 Pack
D0433B-4	125mL Natural Wide Mouth Laboratory Grade Bottle, 38mm Closure - Bulk Pack	350 Pack
D0433-8	250mL Natural Wide Mouth Laboratory Grade Bottle, 45mm Closure	24 Pack
D0433B-8	250mL Natural Wide Mouth Laboratory Grade Bottle, 45mm Closure - Bulk Pack	180 Pack
D0433-16	500mL Natural Wide Mouth Laboratory Grade Bottle, 55mm Closure	12 Pack
D0433B-16	500mL Natural Wide Mouth Laboratory Grade Bottle, 55mm Closure - Bulk Pack	100 Pack
D0433-32	1000mL Natural Wide Mouth Laboratory Grade Bottle, 65mm Closure	12 Pack
D0433B-32	1000mL Natural Wide Mouth Laboratory Grade Bottle, 65mm Closure - Bulk Pack	55 Pack

High Density Polyethylene (HDPE), Amber - Class 1

Cat. No.	Description	Qty
D0434-4	125mL Amber Wide Mouth Laboratory Grade Bottle, 38mm Closure	24 Pack
D0434B-4	125mL Amber Wide Mouth Laboratory Grade Bottle, 38mm Closure - Bulk Pack	72 Pack
D0434-8	250mL Amber Wide Mouth Laboratory Grade Bottle, 45mm Closure	24 Pack
D0434B-8	250mL Amber Wide Mouth Laboratory Grade Bottle, 45mm Closure - Bulk Pack	72 Pack
D0434-16	500mL Amber Wide Mouth Laboratory Grade Bottle, 55mm Closure	12 Pack
D0434B-16	500mL Amber Wide Mouth Laboratory Grade Bottle, 55mm Closure - Bulk Pack	48 Pack
D0434-32	1000mL Amber Wide Mouth Laboratory Grade Bottle, 65mm Closure	12 Pack
D0434B-32	1000mL Amber Wide Mouth Laboratory Grade Bottle, 65mm Closure - Bulk Pack	24 Pack

Laboratory Grade Plastic Bottles

Standard

Narrow Mouth Plastic Laboratory Grade Bottles with Closures

High Density Polyethylene (HDPE), Natural - Class 1

Cat. No.	Description	Qty
D0437-4	125mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0437B-4	125mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	350 Pack
D0437-8	250mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0437B-8	250mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	180 Pack
D0437-16	500mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0437B-16	500mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	100 Pack
D0437-32	1000mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0437B-32	1000mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	55 Pack

High Density Polyethylene (HDPE), Amber - Class 1

Cat. No.	Description	Qty
D0438-4	125mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0438B-4	125mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	72 Pack
D0438-8	250mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0438B-8	250mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	72 Pack
D0438-16	500mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0438B-16	500mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	48 Pack
D0438-32	1000mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0438B-32	1000mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	24 Pack

