

Narrow Mouth Septum Bottles



- Preassembled with PTFE/Silicone-lined open top black polypropylene closures
- Available in standard, precleaned, or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Open Top Black Polypropylene (PP) Closures and Septum

Clear Glass Narrow Mouth Septum Bottles - Class 1

Cat. No.	Description	Qty
9-160	1 oz, 30mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-161	2 oz, 60mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-162	4 oz, 125mL Clear Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-163	8 oz, 250mL Clear Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Amber Glass Narrow Mouth Septum Bottles - Class 1

Cat. No.	Description	Qty
9A-160	1 oz, 30mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-161	2 oz, 60mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-162	4 oz, 125mL Amber Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-163	8 oz, 250mL Amber Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Narrow Mouth Septum Bottles

Precleaned

Assembled with Open Top Black Polypropylene (PP) Closures and Septum

Clear Glass Narrow Mouth Septum Bottles - Class 2

Cat. No.	Description	Qty
9-160-2	1 oz, 30mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-161-2	2 oz, 60mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-162-2	4 oz, 125mL Clear Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-163-2	8 oz, 250mL Clear Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Amber Glass Narrow Mouth Septum Bottles - Class 2

Cat. No.	Description	Qty
9A-160-2	1 oz, 30mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-161-2	2 oz, 60mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-162-2	4 oz, 125mL Amber Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-163-2	8 oz, 250mL Amber Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Precleaned and Certified

Assembled with Open Top Black Polypropylene (PP) Closures and Septum

Clear Glass Narrow Mouth Septum Bottles - Class 3

Cat. No.	Description	Qty
9-160-3	1 oz, 30mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-161-3	2 oz, 60mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-162-3	4 oz, 125mL Clear Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-163-3	8 oz, 250mL Clear Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Amber Glass Narrow Mouth Septum Bottles - Class 3

Cat. No.	Description	Qty
9A-160-3	1 oz, 30mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-161-3	2 oz, 60mL Amber Septum Bottle, 20-40mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-162-3	4 oz, 125mL Amber Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-163-3	8 oz, 250mL Amber Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

