

## Vista Vials™, 12x32mm, Patented

### Vials



The Vista Vial™ Series is a two part design solution combining a glass base with a polyethylene top – just snap the two pieces into place. For the glass base, choose from clear or amber

Type I borosilicate glass. The tops are made of polyethylene and are available in a variety of neck finishes – 8-425mm thread, 9mm thread, 11mm crimp. The Vista Vial™ is easier to fill – the wide opening of the glass bottom facilitates direct pipetting and is ideal for viscous fluids such as detergents and oils. Depending on the sample, disposal of vials is easier as the plastic can be incinerated intact without the need for closure removal.



2.0mL Vista Vial™,  
9mm Thread



2.0mL Vista Vial™,  
11mm Crimp

#### Glass Base (1.0mL volume)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
4115-1220	Clear Glass Vista Vial™ Base
4115-1220A	Amber Glass Vista Vial™ Base

#### Polyethylene Top

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
31508PE-12	8-425mm Thread Polyethylene Vista Vial™ Top
31509PE-12	9mm Thread Polyethylene Vista Vial™ Top
31511PE-12	11mm Snap Seal™ Polyethylene Vista Vial™ Top
31508PE-12A	8-425mm Thread Amber Polyethylene Vista Vial™ Top
31509PE-12A	9mm Thread Amber Polyethylene Vista Vial™ Top
31511PE-12A	11mm Snap Seal™ Amber Polyethylene Vista Vial™ Top

### Closures



Available in a wide assortment of colors. Designed for use with all autosamplers.

#### Screw Thread Closures (8-425mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
5310-08	Black, Polypropylene Large Open Hole Cap for Shimadzu
5310L-08	Black, Polypropylene Large Open Hole Cap
5320-08	Solid Top, Black Polypropylene Unlined
5330-08	Top Seal™ Closure, Clear Polypropylene, molded 10mil septa
5360-08	Solid Top, Black Polypropylene, PTFE/F217 Lined

## Vista Vials™, 12x32mm, Patented

### Closures



Refer to Septa Selection Guide to choose 8mm septa best suited for your application.

#### Septa for Screw Thread Closures (8-425mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
601010-08	8mm x 0.010" White PTFE Septa
603570-08	8mm x 0.035" Red PTFE/Silicone Septa with Slit
604040-08	8mm x 0.045" PTFE/Red Rubber Septa
604040-08F	8mm x 0.045" PTFE/Butyl Rubber Septa
604060-08	8mm x 0.040" Red PTFE/Silicone/Red PTFE Septa
604550-08	8mm x 0.045" Red PTFE/Silicone Septa
606550-08	8mm x 0.065" Red PTFE/Silicone Septa



Choose the preassembled caps with septa best suited for your application. The patented Top Hat™ septa and closure was designed to provide a positive physical lock to assure that

preassembled components arrive intact.

#### Preassembled Screw Thread Closure and Septa (8-425mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
804040-08	Black Closure, PTFE/Red Rubber Septa, 0.040"
804060-08	Black Closure, PTFE/Silicone/PTFE Septa, 0.040"
806550-08	Black Closure, Red PTFE/Silicone Septa, 0.065"
806550TH-08	Black Top Hat™ PTFE/Silicone Septa/Closure (Patented)



Easy-on, easy-off convenience with a twist of the cap that produces an evaporation-proof seal. Larger opening provides increased target area

for needle penetration. Available in a variety of colors for easy sample identification. Choose a septa best suited for your application.

#### Screw Thread Closures (9mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
5310F-09RB	9mm R.A.M.™ Royal Blue, Polypropylene Open Hole Cap
5310-09FRB	9mm R.A.M.™ Ribbed Royal Blue, Polypropylene Open Hole Cap
5320-09	9mm R.A.M.™ Solid Top, Black Polypropylene Unlined Cap
5330-09	9mm R.A.M.™ Top Seal™ Closure, Clear Polypropylene, molded 10mil septa
5360-09	9mm R.A.M.™ Solid Top, Black Polypropylene PTFE/F217 Lined Cap

#### Septa for R.A.M. Screw Thread Closures (9mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
601010-09	9mm x 0.010" White PTFE Septa
604040-09	9mm x 0.040" PTFE/Red Rubber Septa
604040-09F	9mm x 0.040" PTFE/Butyl Rubber Septa
604050-09	9mm x 0.040" PTFE/Silicone Septa
604060-09	9mm x 0.040" Red PTFE/Silicone/Red PTFE Septa
604070-09	9mm x 0.040" Red PTFE/Silicone w/Slit Septa

**J.G. Finneran Associates, Inc.** — Visit [www.JGFinneran.com](http://www.JGFinneran.com) • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • [sales@novatech-usa.com](mailto:sales@novatech-usa.com) • [www.novatech-usa.com](http://www.novatech-usa.com)

# Vista Vials™, 12x32mm, Patented

## Closures



9mm R.A.M.™ Cap

Preassembled closures and septa come ready to use. Designed to work in Agilent or other robotic arm autosamplers. Available in a variety of colors.



9mm R.A.M.™ Ribbed Cap

### Preassembled R.A.M.™ Closure and Septa (9mm)

**Shelf Pack - 100 pieces; Case Pack - 100 pieces**

Cat. No.	Description
5394F-09RB	9mm R.A.M.™ Cap, Royal Blue, PTFE/Butyl Rubber
5395F-09RB	9mm R.A.M.™ Cap, Royal Blue, PTFE/Silicone
53951F-09RB	9mm R.A.M.™ Cap, Royal Blue, Fixed PTFE/Silicone
5396F-09RB	9mm R.A.M.™ Cap, Royal Blue, PTFE/Silicone/PTFE
5397F-09RB	9mm R.A.M.™ Cap, Royal Blue, PTFE/Silicone w/Slit
53971F-09RB	9mm R.A.M.™ Cap, Royal Blue, Fixed PTFE/Silicone w/Slit
5394-09FRB	9mm R.A.M.™ Ribbed Cap, Royal Blue, PTFE/Butyl Rubber
5395-09FRB	9mm R.A.M.™ Ribbed Cap, Royal Blue, PTFE/Silicone
53951-09FRB	9mm R.A.M.™ Ribbed Cap, Royal, Fixed PTFE/Silicone
5396-09FRB	9mm R.A.M.™ Ribbed Cap, Royal Blue, PTFE/Silicone/PTFE
5397-09FRB	9mm R.A.M.™ Ribbed Cap, Royal Blue, PTFE/Silicone w/Slit
53971-09FRB	9mm R.A.M.™ Ribbed Cap, Royal Blue, Fixed PTFE/Silicone w/Slit

### Poly Crimp™ Seals (11mm) - Patented

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
5500-11	11mm Clear Poly Crimp™ Seal, w/Molded Septum
5500SB-11	11mm Clear Poly Crimp™ Seal, w/Starburst
5510-11	11mm Clear Poly Crimp™ Seal, 10mil PTFE
5540-11	11mm Clear Poly Crimp™ Seal, PTFE/Red Rubber
5550-11	11mm Clear Poly Crimp™ Seal, PTFE/Silicone
5560-11	11mm Clear Poly Crimp™ Seal, PTFE/Silicone/PTFE
5570-11	11mm Clear Poly Crimp™ Seal, PTFE/Silicone w/Slit
5581-11	11mm Clear Poly Crimp™ Seal, PTFE/Silicone w/Starburst



A crimp cap constructed of plastic that fits virtually any 12x32mm sample vial with a crimp style finish. Available in a variety of colors for easy sample identification.

**J.G. Finneran Associates, Inc.** — Visit [www.JGFinneran.com](http://www.JGFinneran.com) • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2008 certified guidelines

(866) 433-6682 • (281) 359-8538 • [sales@novatech-usa.com](mailto:sales@novatech-usa.com) • [www.novatech-usa.com](http://www.novatech-usa.com)

## Vista Vials™, 12x32mm, Patented

### Closures



Caps may be snapped on and off by hand which eliminates the need for metal crimping and decapping tools, thus reducing

potential hazards due to sharp aluminum metal edges. Large opening top provides greater target area for needle penetration. Available in a variety of colors for easy sample identification.

#### Snap Top Caps™ (11mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
5200-11	11mm Clear Snap Cap Seal w/molded Septum
5200SB-11	11mm Clear Snap Cap Seal w/Starburst
5210-11	11mm Clear Snap Cap, 10mil PTFE
5240-11	11mm Clear Snap Cap, PTFE/Butyl Rubber
5250-11	11mm Clear Snap Cap, PTFE/Silicone
5260-11	11mm Clear Snap Cap, PTFE/Silicone/PTFE
5270-11	11mm Clear Snap Cap, PTFE/Silicone w/Slit
5275-11GY	11mm Gray Snap Cap, PTFE/Silicone w/Cross
5280-11	11mm Snap Cap, PTFE/Silicone/PTFE w/Starburst
5281-11	11mm Clear Snap Cap, PTFE/Silicone w/Starburst
5282-11	11mm Clear Snap Cap, 10mil PTFE w/Starburst



Available in a variety of colors for easy sample identification.

#### Aluminum Seals (11mm)

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

Cat. No.	Description
5110-11	11mm Silver Seal, 10mil PTFE
5120-11	11mm Silver Seal, 10mil Polypropylene
5140-11	11mm Silver Seal, Clear PTFE/Red Rubber
5140-11X	11mm Silver Seal, PTFE/Natural Red Rubber
5150-11	11mm Silver Seal, PTFE/Silicone
5160-11	11mm Silver Seal, PTFE/Silicone/PTFE
5170-11	11mm Silver Seal, PTFE/Silicone w/Slit
5190-11	11mm Silver Seal, Black Viton®

### Accessories




Snap Racks™ are constructed of durable glass reinforced polypropylene, with a center well to conveniently hold caps, seals, septa, etc. Racks snap together securely to double your rack capacity.

#### Snap Racks™

**Case - 5 pieces**

Cat. No.	Description
9700-12	50 Position White Rack for 12mm Vials
9700-12B	50 Position Blue Rack for 12mm Vials



[sales@novatech-usa.com](mailto:sales@novatech-usa.com)  
[www.novatech-usa.com](http://www.novatech-usa.com)  
Tel: (866) 433-6682 Fax: (866) 433-6684  
Tel: (281) 359-8538 Fax: (281) 359-0084