

Wide Mouth Septum Jars



-Preassembled with ultrasonically welded PTFE/Silicone liner

-Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Prcleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Open Top White Polypropylene Closures and Septum

Tall Straight Sided Wide Mouth Septum Jars - Class 1

Cat. No.	Description	Qty
9-141	4 oz, 125mL Clear Tall Septum Jar, 48-400mm Closure, PTFE/Silicone Liner	24 Pack
9-142	8 oz, 250mL Clear Tall Septum Jar, 58-400mm Closure, PTFE/Silicone Liner	12 Pack

Short Straight Sided Wide Mouth Septum Jars - Class 1

Cat. No.	Description	Qty
9-140	2 oz, 60mL Clear Short Septum Jar, 53-400mm Closure, PTFE/Silicone Liner	24 Pack

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
9-150	500cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack
9-151	950cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack

Wide Mouth Septum Jars

Precleaned

Assembled with Open Top White Polypropylene Closures and Septum

Tall Straight Sided Wide Mouth Septum Jars - Class 2

Cat. No.	Description	Qty
9-141-2	4 oz, 125mL Clear Tall Septum Jar, 48-400mm Closure, PTFE/Silicone Liner	24 Pack
9-142-2	8 oz, 250mL Clear Tall Septum Jar, 58-400mm Closure, PTFE/Silicone Liner	12 Pack

Short Straight Sided Wide Mouth Septum Jars - Class 2

Cat. No.	Description	Qty
9-140-2	2 oz, 60mL Clear Short Septum Jar, 53-400mm Closure, PTFE/Silicone Liner	24 Pack

Amber Glass Wide Mouth Packers - Class 2

Cat. No.	Description	Qty
9-150-2	500cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack
9-151-2	950cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack

Precleaned and Certified

Assembled with Open Top White Polypropylene Closures and Septum

Tall Straight Sided Wide Mouth Septum Jars - Class 3

Cat. No.	Description	Qty
9-141-3	4 oz, 125mL Clear Tall Septum Jar, 48-400mm Closure, PTFE/Silicone Liner	24 Pack
9-142-3	8 oz, 250mL Clear Tall Septum Jar, 58-400mm Closure, PTFE/Silicone Liner	12 Pack

Short Straight Sided Wide Mouth Septum Jars - Class 3

Cat. No.	Description	Qty
9-140-3	2 oz, 60mL Clear Short Septum Jar, 53-400mm Closure, PTFE/Silicone Liner	24 Pack

Amber Glass Wide Mouth Packers - Class 3

Cat. No.	Description	Qty
9-150-3	500cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack
9-151-3	950cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack

