

LabStrong Replacement Cartridge Pack for Nanopure® Diamond™ Water System

LabStrong Cartridge Packs can be used in the Nanopure® Diamond™ water purifications systems to efficiently produce ASTM Type I water without any leaks!








Product Features

- Purity up to 18.2 megohms-cm
- Leak free technology
- Uses two types of activated carbon to absorb organic compounds and chlorine, and the highest purity semiconductor-grade ion-exchange resins to remove final ionic contaminants. Systems fed with pre-treated water.
- Variety of four packs to choose from
- Easy Installation with one connection
- Tested to meet ASTM Type I Standards
- Meets strictest requirements for TOC and Resistivity

Plus! All Packs include one 0.2µm final filter

**100% Leak
Free Tested**

 (Water) Every Drop Matters... 

LabStrong #	Barnstead #	Description	Where Used	Price
D50280LS	D50280	Cartridge Pack Organic Free R/O and Distilled Feed 	Ultra-Low Organics Type I Water Reverse Osmosis or Distilled Water Feed	\$797.00
D50281LS	D50281	Cartridge Pack Organic Free Deionized Water Feed 	Ultra-Low Organics Type I Water Deionized Water Feed	\$797.00
D50282LS	D50282	Cartridge Pack Type I R/O and Distilled Feed 	Low Organics Type I Water Reverse Osmosis or Distilled Water Feed	\$797.00
D50283LS	D50283	Cartridge Pack Type I Deionized Water Feed 	Low Organics Type I Water Deionized Water Feed	\$797.00
CMX25LS	CMX25	Cleaning Syringe 	EASypure II (LF, UV & UV/UF), NANOpure Infinity & NANOpure Diamond	\$36.00
FL0016-1	D3750	Sterile 0.2 Micron Final Filter with Endbell	B-Pure, E-Pure, all NANOpure systems (including Diamond), AccuDispense, EASypure, all EASypure II systems, EASypure RoDi & Diamond TII	\$99.00
LMP13	LMX13	185/254 nm UV Lamp	NANOpure Quad, Infinity, EASypure, EASypure II (RF/UV, UV & UV/UF), EASypure RoDi, older NANOpure Diamond & Diamond TII	\$138.00
LMP01	LMX13	185/254 nm UV Lamp	Diamond TII & NANOpure Diamond with new style ballast	\$138.00
PU1190X1LS	PU1190X1	Diaphragm Pump	NANOpure Infinity Retrofit Kit (PU896X3) & NANOpure Diamond	\$399.00

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084

MADE IN USA WITH GLOBAL MATERIALS

 **LabStrong**TM
Strength Through Innovation

LabStrong #	Barnstead #	Description	Where Used	Price
D0803LS	D0803	High Capacity (Two-Bed)	B-Pure, Bantam, Cooling Loops, E-Pure, NANOpure, NANOpure II & NANOpure Quad	\$132.00
D0809LS	D0809	Ultrapure (Mixed-Bed)	B-Pure, Bantam, Cooling Loops, E-Pure, NANOpure, NANOpure II & NANOpure Quad	\$136.00
D0813LS	D0813	Organic Removal (Carbon)	B-Pure, Bantam & Cooling Loops	\$111.00
D0815LS	D0815	Cation Removal (Cation) NEW!	B-Pure, Bantam & Cooling Loops	\$107.00
D0820LS	D0820	Organic Free	NANOpure II & older NANOpure	\$164.00
D0832LS	D0832	Pretreatment Cartridge (Mixed-Bed/Carbon)	B-Pure & Still Pretreatment	\$120.00
D0835LS	D0835	Combination Cartridge (Colloid/Organic)	E-Pure, older NANOpure, NANOpure II & NANOpure Quad	\$131.00
D0836LS	D0836	Macropure (Macroreticular/Carbon)	E-Pure, older NANOpure, NANOpure II & NANOpure Quad	\$174.00
D8901LS	D8901	HN (3/8") Dyed High Capacity (Two-Bed)	Hose Nipple, non-pressure application	\$130.00
D8902LS	D8902	HN (3/8") Dyed Ultrapure (Mixed-Bed)	Hose Nipple, non-pressure application	\$130.00
D8904LS	D8904	HN (3/8") Organic Removal - No Dye (Carbon)	Hose Nipple, non-pressure application	\$106.00
D8905LS	D8905	HN (3/8") Dyed Cation Removal NEW!	Hose Nipple, non-pressure application	\$105.00
D8911LS	D8911	HN (3/8") Ultrapure - No Dye (Mixed-Bed)	Hose Nipple, non-pressure application	\$131.00
D8921LS	D8921	HN (3/8") Dyed Pretreatment (Mixed-Bed/Carbon) NEW!	Hose Nipple, non-pressure application	\$131.00
D8922LS	D8922	HN (3/8") Dyed Organic & Ultrapure Removal (Mixed-Bed/Carbon)	Hose Nipple, non-pressure application	\$140.00
D8950LS	D8950	HN (3/8") 2/3 Size Dyed High Capacity (Two-Bed) NEW!	Hose Nipple, non-pressure application	\$112.00
D8951LS	D8951	HN (3/8") 2/3 Size Dyed Organic Removal (Mixed-Bed/Carbon) NEW!	Hose Nipple, non-pressure application	\$120.00
D5021LS	D5021	Organic Free	E-Pure & NANOpure Quad	\$171.00
D5027LS	D5027	Ultrapure (Mixed-Bed)	E-Pure, NANOpure Quad & ROPure ST	\$134.00
D5028LS	D5028	E-Pure Cartridge Kit (Tap Feed/4 Mod Standard)	E-Pure	\$557.00
D5023LS	D5023	E-Pure Cartridge Kit (Tap Feed/4 Mod Organic Free)	E-Pure	\$642.00
D50227LS	D50227	E-Pure Cartridge Kit (Pretreated Feed/4 Mod Standard)	E-Pure	\$559.00
D50228LS	D50228	E-Pure Cartridge Kit (Pretreated Feed/4 Mod Organic Free)	E-Pure	\$642.00
D5029LS	D5029	E-Pure Cartridge Kit (3 Mod Standard)	E-Pure	\$429.00
D5022LS	D5022	E-Pure Cartridge Kit (3 Mod Organic Free)	E-Pure	\$615.00
D4801LS	D4801	NANOpure Quad Style Cartridge Kit (Tap Feed/Standard)	NANOpure Quad	\$589.00
D4802LS	D4802	NANOpure Quad Style Cartridge Kit (Tap Feed/Organic Free)	NANOpure Quad	\$731.00
D5025LS	D5025	NANOpure Quad Style Cartridge Kit (Pretreated Feed/Standard)	NANOpure Quad	\$569.00
D5026LS	D5026	NANOpure Quad Style Cartridge Kit (Pretreated Feed/Organic Free)	NANOpure Quad	\$772.00
D3804LS	D3804	NANOpure II Cartridge Kit (4 Mod Standard)	NANOpure II & Older NANOpure	\$629.00
D3805LS	D3805	NANOpure II Cartridge Kit (3 Mod Standard)	NANOpure II & Older NANOpure	\$489.00
D3806LS	D3806	NANOpure II Cartridge Kit (3 Mod Organic Free)	NANOpure II & Older NANOpure	\$549.00
D3807LS	D3807	NANOpure II Cartridge Kit (4 Mod Organic Free)	NANOpure II & Older NANOpure	\$699.00
D400377LS	D400377	Mega-Pure High Purity Cartridge NEW!	Mega-Pure D1, D2 & MP12-A Stills	\$249.00
D400499LS	D400499	Mega-Pure High Capacity Cartridge NEW!	Mega-Pure D1, D2 & MP12-A Stills	\$259.00
D440265LS	D440265	Mega-Pure Organic and Chlorine Removal Cartridge NEW!	Mega-Pure D2 & MP12-A Stills	\$232.00
06808LS	06808 & 06808-DB	NANOpure II Canister O-Ring NEW!	NANOpure II	\$6.00
25001LS	25001 & 25001-DB	Air Filter	Water Storage Reservoirs (Cabinetized Stills, Classic Stills, Diamond TII, etc.)	\$65.00
410671LS	410671	Mega-Pure Gasket 4 Pack NEW!	Mega-Pure Single (D1) & Dual (D2) Deionizer Systems & Mega-Pure MP-12A Stills	\$69.00
740880LS	740880	Mega-Pure MP-1 Heating Element NEW!	Mega-Pure Stills (MP-1)	\$595.00
740883LS	740883 & 740895	Mega-Pure Heating Element	Mega-Pure Stills (MP-3A, MP-6A, MP-11A & MP-12A)	\$662.00
D8900LS	D8900	Wall Mount Bracket	Hose Nipple Cartridges	\$42.00
FL1265X1LS	FL1265X1	Encapsulated RO Membrane	Diamond RO & Diamond TII	\$472.00
GSX27LS	GSX27	Cartridge O-Ring	B-Pure, E-Pure & NANOpure Quad	\$4.00
GSX28LS	GSX28	Canister O-Ring	B-Pure, E-Pure & NANOpure Quad	\$6.00

*Barnstead, E-Pure, NANOpure, B-Pure, Bantam, EASYPure, Mega-Pure, Thermo Scientific & Diamond are trademarks of Thermo Fisher Scientific. LabStrong cartridges and parts are not affiliated with, nor supplied or sponsored by Thermo Fisher Scientific.



Need other replacement parts for your unit?
 Visit our website for more information on products and replacement parts.

LT0016-7 rev 2 | 2016

(866) 433-6682 • (281) 359-8538 • sales@novatech-usa.com • www.novatech-usa.com