


**LEVGO, INC.**

**We make smartTools for the Lab.™**

## Eco-smartFunnel™ Biodegradable Laboratory Funnel



1/4" Opening

5 inches



The mark of responsible forestry



Flat folded funnels store easily and remain clean until used.  
25 Flat folded funnels/Pack



Mini : 2 3/8" Top Opening  
5" Tall x 1/4" in<sup>2</sup> Stem Opening  
6.3 Oz, 186 mL Est. Volume

### The Eco-IN and Eco-OUT Alternative to Single Use Plastic Funnels

The LevGo Eco-smartFunnel is Eco-In because it is made of FSC® Certified paper and Eco-Out because it is biodegradable in the standard municipal waster stream.

- Innovative "backrest" allows hands-free use when pouring viscous liquids
- Fits a wide variety of openings
- Maximum flexibility - cut stem to any desired opening size or shape
- Apply a minimum of pressure and the funnel will stay in position
- Use with virtually any liquid or powder\*

\*not recommended for hot liquids

[www.LevGo.com](http://www.LevGo.com)

Style	Features	Dimensions	Pack Size	Cat. No.
White Glossy Food Grade Card Stock	Adjustable Stem, Flat Folded, Biodegradable	Top: 2.4 in <sup>2</sup> , 6 cm <sup>2</sup> Stem: 0.25 in <sup>2</sup> , 0.64 cm <sup>2</sup>	25/PK	19630