

LEVGO, INC.

We make smartTools for the Lab.™

smartBoats™ Biodegradable Paper Weigh Boats

LevGo's paper weighing boats are a first in Eco Friendly Lab Consumables

They are both Eco-In and Eco-Out.



The mark of
responsible forestry

Eco-In: Because they are FSC® Certified. The papers used to manufacture our boats are made from responsibly forested wood products.

Eco-Out: Because they are naturally biodegradable in the standard municipal waste stream.

- smartBoats control contents like a polystyrene weigh boat and are easier to handle than standard flat weighing paper.
- Made of grease proof paper that meets FDA 21 CFR 176.170, components of paper in contact with aqueous and fatty foods, and with 21 CFR 176.180 for contact with dry foods.
- Compatible with heat applications with a temperature limitation of 210 °C.
- Small brown & white and large brown & white are made in Italy. Medium white and medium brown are made in the USA.



Medium Boats



www.LevGo.com

Style	Features	Dimensions	Pack Size
Small Brown	Biodegradable, FSC® Certified	1⅜ in² x ⅝ in, 3.5 cm² x 16 cm	500/PK
Medium Brown	Biodegradable, FSC® Certified	1⅞ in² x ⅞ in, 4.8 cm² x 14 cm	500/PK
Large Brown	Biodegradable, FSC® Certified	3¾ in² x 1 in, 9.5 cm² x 2.5 cm	NEW! 250/PK
Small White	Biodegradable, FSC® Certified	1⅜ in² x ⅝ in, 3.5 cm² x 16 cm	500/PK
Medium White	Biodegradable, FSC® Certified	1⅞ in² x ⅞ in, 4.8 cm² x 14 cm	500/PK
Large White	Biodegradable, FSC® Certified	3¾ in² x 1 in, 9.5 cm² x 2.5 cm	250/PK