


LEVGO, INC.

We make smartTools for the Lab.™

smartSlicer™ Gel Cutter and Laboratory Safty Razor




At just 0.27 mm (0.011") the VWR Gel Cutter is as thin as a single edge industrial razor blade but much safer

- the LevGo smartSlicer is unlikely to cut you!

- The LevGo smartSlicer is able to cleanly slice gels right on the light Box without scratching the surface of the box.*
- The LevGo smartSlicer can also be used to chop, mince or slice soft tissue samples.
- The LevGo smartSlicer is also adjustable - simply cut it with standard scissors to the width you require.
- Available in two sizes: Regular, sized the same as an industrial single edge razor and Large, 30% larger for bigger jobs.
- LevGo smartSlicers are made in the USA of food grade PET and packed 100 per pack.
- Easy to sterilize with alchohol or treat with DNase or RNase cleaning agents.

*Best for 1% to 2% gels.



Lab Safety Managers love the LevGo smartSlicer!

www.LevGo.com

Style	Features	Dimensions	Pack Size	Cat. No.
Small (same as regular razor)	Handy grip area, cut to size	¾" x 1.5" x 0.01" , 19mm L x 38mm H , 0.27mm thick	100/PK	30202
Large	Handy grip area, cut to size	1" x 2" x 0.01" , 25mm L x 50mm H , 0.27mm thick	100/PK	30203