



We make smartTools for the Lab.™


# smartPicks™ Biodegradable Inoculation Picks

Based on a long tradition of using Toothpicks in the Lab.


## Pick



## Drop



## Close



smartPicks Mini Wooden Inoculation Picks

A faster way to inoculate.

- Autoclave
- **Pick** your colony
- **Drop** the Mini Pick into your liquid
- **Close** the cap
- No pick removal step.

Available in Semi (blunt) or needle point.

## Mini



1.28" Tall  
0.1" diameter

- Made in the USA in an ISO 9001 : 2008, ISO 13485 : 2003 & FDA registered facility
- 100% USA grown White Birch
- Non-Sterile, smartPicks meet FDA requirements for food
- **Autoclavable**
- Environmentally Friendly - Naturally Biodegradable



## Standard



2.5" Tall  
0.1" diameter

## Large



6" Tall  
0.125" diameter

Available in 3 convenient sizes, 1.28", 2.5" and 6" to accommodate various transfer, inoculating and streaking needs.

[www.LevGo.com](http://www.LevGo.com)

Description	Features	Dimensions (L)	Pack Size	Cat No.
Mini	Semi Point, Fits in 1.5ml tube	1.28", 3.25 cm	2000/PK	18128SP
Mini	Needle Point, Fits in 1.5ml tube	1.28", 3.25 cm	2000/PK	18128NP
Standard	Semi Point	2.5", 6.4 cm	1600/PK	18250SP
Standard	Needle Point	2.5", 6.4 cm	1600/PK	18250NP
Large	Semi Point	6", 15.2 cm	500/PK	18600SP
Large	Needle Point	6", 15.2 cm	500/PK	18600NP