

BluBlot™ HS AUTORADIOGRAPHY FILM

novatech
INTERNATIONAL

sales@novatech-usa.com


www.novatech-usa.com

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084



- **High Contrast**
- **Reduced Background**
- **Ideal for Chemiluminescence**
- **Publication Quality**
- **Suitable for ¹⁴C, ³²P, ³³P, ¹²⁵I, ³⁵S**
- **Double Emulsion**

BluBlot HS is a high quality, blue sensitive film that is ideal for a wide range of electrophoretic and autoradiography applications. It can be used with calcium tungstate or blue-emitting rare-earth screens.

BluBlot HS exhibits excellent sensitivity and contrast characteristics and is coated on both sides for improved consistency and overall performance. It reliably produces sharp, bright images that are easily visualized and suitable for publication.

Especially useful with chemiluminescence, BluBlot HS works well with ¹⁴C, ³²P, ³³P, ¹²⁵I, and ³⁵S.

Compatible with automatic processors and manual development methods, BluBlot HS compares directly with films costing up to 3 times as much.

With a shelf life of 2 years from date of production, BluBlot HS is packaged in lined, light-safe dispenser boxes of 100 sheets per box.

MTCBio

2600A Main Street Extension

Sayreville, NJ 08872

Ph: 908-279-0280 www.mtcbio.com

Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Quantity</u>
A8803	BluBlot HS Autoradiography film, 5x7 (127x178mm)	100 Sheets
A8805	BluBlot HS Autoradiography film, 8x10 (203x254mm)	100 Sheets