

# MTC Bio™ INOCULATING LOOPS & NEEDLES



sales@novatech-usa.com  
www.novatech-usa.com  
Tel: (866) 433-6682 Fax: (866) 433-6684  
Tel: (281) 359-8538 Fax: (281) 359-0084



- **Sterile & disposable... no flaming**
- **Uniform streaking & inoculating**
- **Hexagonal grip - improves control**
- **Smooth loop surfaces - no gouging**
- **Nontoxic, biologically inert polypropylene**

*These sterile disposable plastic loops and needles are designed for use in quantitative procedures such as sampling, urine counts and serial dilutions, as well as for bacterial inoculation. Because they are presterilized and designed for single use, there is no risk of cross contamination and no need for flaming. They are flexible enough for streaking, yet rigid enough for subsurface inoculation.*

*The unique hexagonal shape of the shaft makes it easier to grip, orientate and guide through a smooth streaking process. The loop surfaces are extremely smooth with no rough edges and no flashing or burrs on the loop head. This results in the elimination of cutting or gouging of the agar surface.*

*Sterile, disposable and verifiably accurate, these loops provide an efficient, cost-effective and time-saving alternative to traditional metallic loops.*



2600A Main Street Extension  
Sayreville, NJ 08872  
Ph: 908-279-0280 www.mtcbio.com

## Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Length</u>	<u>Packaging</u>	<u>Quantity</u>
M4015	Inoculating loops, 1ul	200mm	25 per sterile zip bag	1,000
M4025	Inoculating loops, 10ul	200mm	25 per sterile zip bag	1,000
M4030	Inoculating needles	200mm	25 per sterile zip bag	1,000