

LABeler™ Handheld Laboratory Printer



- **Lightweight, handheld operation**
- **Variety of label tapes for many tube sizes**
- **Built in cutter - no need for scissors**
- **Prints a variety of barcodes**
- **“Ultra-Low” tape for liquid nitrogen storage**

Identifying and organizing increasingly large numbers of samples has become an important objective for today's busy research laboratory. The process of writing on sample tubes & containers can be messy and unreliable. It can also be difficult to find the time to access a computer and design a layout for expensive ink jet labels. The new LABeler™ Lab Printer provides the convenience of a lab marking pen with the benefits of a large desktop printer. Simply type in the desired information using the 50+ button key pad and press the print key...just press the cut lever and your label is instantly ready for application.


With over 12 available sizes, colors & material combinations, the LABeler is ideal for label identification of any laboratory vessel for nearly any laboratory application. Standard label tapes can be autoclaved, boiled or even frozen to temperatures as low as -50°C. Special Ultra-Low label tapes are also available for storage in -196°C environments or for liquid nitrogen applications



Available with 9V power supply
(also runs off 6 AA batteries)



Replaceable tape cassettes
(26ft / 792cm length)



2600A Main Street Extension
Sayreville, NJ 08872
Ph: 908-279-0280 • www.mtcbio.com

Ordering Information

Cat. No.

L9010
L9010-A
L9010-PS
L9010-6WK
L9010-6CK
L9010-12WK
L9010-12CK
L9010-12YK
L9010-12RW
L9010-12BW
L9010-12GW
L9010-24WK
L9010-24CK
L9010-12ULT
L9010-24ULT

Description

LABeler™ Handheld lab printer, with 9V power supply (U.S. Plug)*
LABeler™ Handheld lab printer, without power supply
Power supply (9V) for LABeler™ Lab Printer (U.S. Plug)*
26' Cassette of 6mm lab tape, white with black print
26' Cassette of 6mm lab tape, clear with black print
26' Cassette of 12mm lab tape, white with black print
26' Cassette of 12mm lab tape, clear with black print
26' Cassette of 12mm lab tape, yellow with black print
26' Cassette of 12mm lab tape, red with white print
26' Cassette of 12mm lab tape, blue with white print
26' Cassette of 12mm lab tape, green with white print
26' Cassette of 24mm lab tape, white with black print
26' Cassette of 24mm lab tape, clear with black print
20' Cassette of 12mm ULT (ultra-low temp) lab tape, white with black print
20' Cassette of 24mm ULT (ultra-low temp) lab tape, white with black print

* For EU plug, please add "-E" to the end of the item number