

ProPette™ MPC

MANUAL PIPETTE CONTROLLER



- **For glass or plastic pipettes**
- **Accommodates pipettes up to 100ml**
- **Easy, single lever control**
- **Autoclavable lower body**
- **Fingertip speed control**

ProPette™ MPC is the precise and easy to use manual pipette controller that's always ready to use. After attaching a pipette, one squeeze of its large rubber filling bulb is all it takes to start pipetting - there are no wires, no tubing and no battery to recharge! Begin aspirating by touching the EZ-Flow™ lever which is designed to precisely control the filling speed, based on how gently or firmly you press the lever.

A built-in replaceable filter ensures optimum flow rate and protects the controller from liquid penetration and over-fill. A separate button provides a positive blow-out function when required.

Made to accept pipettes from 0.1 to 100ml, the silicone "chuck" adapter is flexible, completely autoclavable, and will accommodate all types of serological, "blow-out", "to deliver," glass and plastic pipettes with a soft, precise, air-tight seal.

MTC Bio Corp.
1 Industrial Way Westfield, NJ 07090
2600A Main Street Extension
Easton, NJ 07724
Sayreville, NJ 08872
732-279-7777 www.mtcbio.com
Ph: 908-279-0280

Ordering Information

<u>Cat. No.</u>	<u>Description</u>
P6420	ProPette MPC, Manual Pipette Controller with filter
P6080-FI	Replacement filter, 0.45um, package of 5