

# MTC Bio<sup>TM</sup> STORAGE RACKS & BOXES

- Polypropylene, withstand -80°C
- 6 types - for 0.2ml to 50ml
- Compact & space efficient
- Autoclavable
- Numbered tube wells

## Test Tube Racks

MTC Bio's benchtop tube racks are made of 100% polypropylene and are available in three popular configurations. No. R1030 is a 4-way rack that can be used to hold 4x 50ml, 12x 15ml or 32x microtubes (0.5, 1.5 or 2.0ml). They can be snap-linked together to form as large a rack as is needed. No. R1070 is a unique rack that holds 30x 15ml and 20x 50ml. No. R1050 is a reversible rack for 96x 1.5ml on one side and 96x 0.5ml on the other.


## Freezer/Storage Boxes

Suitable for refrigerator and freezer storage as low as -80°C, these boxes have translucent lids and numbered cavities. No. R1010 is a PCR workstation rack/box for 96x 0.2ml (plates, strips or individual tubes). No. R1020 is a 10x10 format with an attached hinged lid. No. R1060 is the standard 9x9 format with a friction fit lid.

**MTC Bio**

2600A Main Street Extension  
Sayreville, NJ 08872  
Ph: 908-279-0280 www.mtcbio.com

## Ordering Information



Cat. No.	Description	Dimensions	Color
R1030	4-way rack for 4x50ml, 12x15ml, 32x0.5/1.5ml	6.9x3.7x2.0in	Rainbow Pack
R1070	Combo-rack for 30x15ml and 20x50ml	8.1x6.7x2.4in	Orange
R1050	Reversible rack for 96x1.5ml or 96x0.5ml	9.9x4.4x0.7in	Orange
R1010	Box rack with lid for 96x0.2ml (tubes, strips or plate)	5.0x3.3x1.2in	Orange
R1020	Freezer box, 100x2ml, with hinged lid	5.6x5.8x2.2in	Rainbow Pack
R1060	Freezer box, 81x2ml, with friction fit lid	5.1x5.1x2.0in	Rainbow Pack

Above packs of racks are also available in individual colors as well. Please choose from Red, Blue, Green, Yellow, Orange and Black. Please add the first letter of the desired color to the end of the catalog number. (For example: R1020-R for Red). Add "BK" for Black.