

# SureSeal S™

## 1.5mL MicroTubes

Always Supplied Sterile  
RNase,  
DNase &  
Pyrogen  
Free  
At No Extra Charge

**Syringe Permeable**  
2mm circle is thinner,  
allowing through-cap injection

**Raised locking rim**  
seals into groove in sidewall

**Smooth Rounded Lip**  
prevents "red thumb syndrome"

**Glass-clear side walls**  
maximum transparency  
BioClear™ polypropylene

**Ultra-smooth interior surface**  
Nano-polished molding process

**Advanced BioClear™ polypropylene**  
enhances pellet visualization  
and minimizes optical distortion

**Ring and Groove closure**  
protruding ring on cap locks  
into groove in upper sidewall  
with a distinct "click"

**Frosted marking surface**  
on non-graduated side

**Vivid graduations**  
at 0.1ml and every 0.25ml

**Reinforced Conical Bottom**  
withstands extreme g-forces

- **Sterile comes Standard**  
by electron beam radiation
- **Self-Standing bag**  
with gusseted bottom
- **Stop-Pops™ Supplied**  
10 Stop-Pop cap locks in every bag



Also available in  
these great colors!

MTC Bio's **SureSeal S™** microcentrifuge tubes are designed for a broad range of research applications, including enzyme digests, DNA isolation/purification, centrifugation, incubation and general sample/reagent storage. With a variety of design advances, SureSeal Tubes improve sample visualization, safety, stability and overall ease of use.

**SureSeal S** tubes are manufactured under clean room conditions, using only **BioClear™** co-polymer, an advanced polypropylene that maximizes clarity. **BioClear** has been formulated to be crystal clear, without compromising centrifugal strength (tested up to 28,000xg) or the usable temperature range (-90°C to +121°C).

The oversized tube cap (with frosted marking surface), combined with the unique ring & groove closure, enhances the security of the liquid-tight seal, while minimizing the required locking forces for a more ergonomic opening/closing of the tube caps.

All **SureSeal S** tubes are RNase, DNase and pyrogen free.

**MTC Bio**

2600A Main Street Extension  
Sayreville, NJ 08872

Ph: 908-279-0280 [www.mtcbio.com](http://www.mtcbio.com)

### Ordering Information

Cat. No.	Description	Qty
C2000*	1.5mL microtube, clear, sterile, self-standing bag, 10 Stop-Pops included	500
C2000-AST	1.5mL microtube, sterile, self-standing bag, 10 Stop-Pops, assorted colors	500
C2002	2.0mL microtube, clear, sterile, self-standing bag, 10 Stop-Pops included	500
C2007	0.5mL microtube, clear, sterile, self-standing bag,	500
C2007-AST	0.5mL microtube, sterile, self-standing bag, assorted colors	500

\* Please add (-B) for blue (-R) for red, (-G) for green, (-O) for orange, (-OB) for opaque black, (-P) for purple and (-Y) for yellow