

SureSeal™ Screw-Cap Tubes

With Leak-proof O-ring seal

- **Ultra-clear polypropylene**
- **Graduated at 0.25, 0.50, 0.75, 1.0 and 1.25mL**
- **Leak-proof O-ring seal**
- **RNase, DNase, & Pyrogen free**

SureSeal™ sterile screw-cap microtubes come in a variety of sizes, cap configurations, graduation options, and now light-blocking amber. Made from ultra-clear polypropylene, these screw cap microtubes include assembled caps with silicone O-rings that assure a leakproof seal from -90°C to +121°C. All SureSeal™ screw cap microtubes are sterile, RNase, DNase, & pyrogen free and may be centrifuged at up to 22,000xg.



C2230

The SureSeal™ family includes 1.5mL and 2.0mL tubes with separate or tethered caps and your choice of molded or printed graduations!



C3150-S



C3150-SL



C2230-AM



C3220-SG



C3220-SGL

Cat. No.	Description	Quantity	List Price
C3150-S	SureSeal™ 1.5mL (12.8 x 45.4mm), w/ screw-cap, non-graduated, w/ O-ring, conical, caps assembled, 100 per sterile resealable bag, 10 bags per case	1000/cs	\$178.00
C3150-SL	SureSeal™ 1.5mL (13 x 47.6mm), w/ looped screw-cap, non-graduated, w/ O-ring, conical, caps assembled, 100 per sterile resealable bag, 10 bags per case	1000/cs	\$189.00
C2230	SureSeal™ 2.0mL (12.8 x 46.7mm), w/ screw-cap, molded graduations, w/ O-ring, self-standing, caps assembled, 100 per sterile bag, 10 bags per case	1000/cs	\$178.00
C2230-AM	SureSeal™ 2.0mL (12.8 x 47.4mm), w/ amber screw-cap, molded graduations, w/ O-ring, self-standing, amber caps assembled, 100 per sterile bag, 10 bags per case	1000/cs	\$188.00
C3220-SG	SureSeal™ 2.0mL (13.1 x 46.8mm), w/ screw-cap, printed graduations, w/ O-ring, self-standing, caps assembled, 50 per sterile resealable bag, 20 bags per case	1000/cs	\$229.50
C3220-SGL	SureSeal™ 2.0mL (13 x 49.2mm), w/ looped screw-cap, printed graduations, w/ O-ring, self-standing, caps assembled, 50 per sterile resealable bag, 20 bags per case	1000/cs	\$239.50