

# Weigh Boats

## DISPOSABLE BALANCE DISHES



- **Antistatic and biologically inert**
- **Smooth surface/rounded corners**
- **Recessed corners; precise pouring**
- **Wide, flat bottoms - resist tipping**
- **For liquids, powders or granular samples**
- **White or black for visual contrast**

MTC Bio Weigh Boats are made of antistatic polystyrene and will withstand temperatures to 90°C, as well as exposure to diluted acids, alcohols, bases and aqueous solutions. Their flexibility and smooth interior surfaces with rounded corners allow for easy decanting of powdered samples. Their flat bottoms ensure stability and prevent tipping.

Weigh Boats are available in the traditional square shape as well as in a new diamond shape for easier pouring. Originally designed as weighing containers, they have proven ideal for dispensing & storage and as disposable mixing/quick freezing trays. They are available in black or white for improved sample contrast/visibility and packaged in ultra-compact dispenser packs.



2600A Main Street Extension  
Sayreville, NJ 08872  
Ph: 908-279-0280 www.mtcbio.com

### Ordering Information

#### **Square Weigh Boats**

<u>Cat. No.</u>	<u>Description</u>	<u>Color</u>	<u>Qty.</u>
B6001W	20ml, 45x45mm	White	500
B6002W	100ml, 80x80mm	White	500
B6003W	250ml, 140x140mm	White	500
B6001B	20ml, 45x45mm	Black	500
B6002B	100ml, 80x80mm	Black	500
B6003B	250ml, 140x140mm	Black	500

#### **Diamond Weigh Boats**

<u>Cat. No.</u>	<u>Description</u>	<u>Color</u>	<u>Qty.</u>
B6501W	5ml, 35x55mm	White	500
B6502W	30ml, 65x85mm	White	500
B6503W	100ml, 100x125mm	White	250
B6501B	5ml, 35x55mm	Black	500
B6502B	30ml, 65x85mm	Black	500
B6503B	100ml, 100x125mm	Black	250