

Western Blot Boxes

BLOTTING CONTAINERS



- **A Perfect Box for Every Blot**
- **Reduce Antibody Usage**
- **Hinged Lid Minimizes Evaporation**
- **For Halves, Wholes, Strips, or Minis**



B1206
Multibox™ 6



B1300-8BK

No matter what size your Western blot is, there's a box that fits it perfectly. With boxes that are correctly matched to common blot sizes, you can fit more blots on shakers/rockers while keeping reagent and antibody usage to an absolute minimum. Additionally, all corners are formed at a 90° angle and the bottom surface is uniformly flat, further reducing consumption of costly reagents and antibodies. (For example, popular mini-gels require only about 4 ml.)

At a price that makes them truly disposable, all sizes are engineered from crystal clear polystyrene (except B1200-7BK & B1200-15BK, which are an opaque black that protects light sensitive samples). All sizes are also equipped with a easy-to-open hinged lid that reduces evaporation.

Cat No.	Description	Applications	Minimum Volume
B1200-3	Micro-strip	Various	1 ml
B1200-5	Mini-rectangle	Half mini-protein gels	3 ml
B1200-7	Medium Rectangle	Mini-protein gels, Bio-Rad Protean 3 & similar	4 ml
B1200-9	Long Strip	Protein 3 (cut up)	2 ml
B1200-11	Medium Strip	Various	3 ml
B1200-13	Large	Cambrex PAGEr Gold, Mighty Small, Novex SureLock Mini	6 ml
B1200-15	Square	Owl Penguin P9DS/Puffin	7 ml
B1200-17	Extra Large	Bio-Rad Criterion, Amersham	11 ml
B1200-19	XXL	Bio-Rad Criterion, Amersham	12 ml
B1206	Long Strip (x6)	Protein 3 (multiple)	2 ml x 6
B1300-8BK	Medium Square	Novex Mini Gels	4 ml



2600A Main Street Extension
Sayreville, NJ 08872
Ph: 908-279-0280 www.mtcbio.com

Ordering Information

Cat. No.	Dimensions
B1200-3	2 7/8 x 1 3/16 x 3/4 in (7.3 x 3.0 x 1.9 cm)
B1200-5	2 7/8 x 2 x 1 1/4 in (7.3 x 5.1 x 3.2 cm)
B1200-7	3 1/2 x 2 9/16 x 1 in (8.9 x 6.5 x 2.5 cm)
B1200-7BK	3 1/2 x 2 9/16 x 1 in (8.9 x 6.5 x 2.5 cm)
B1300-8BK	8.6 x 8.6 x 2.8cm (pull off lid)
B1200-9	3 3/4 x 1 3/16 x 5/8 in (9.5 x 3.0 x 1.6 cm)
B1200-11	4 3/8 x 2 1/16 x 1 1/8 in (11.1 x 5.2 x 4.4 cm)

Cat. No.	Dimensions
B1200-13	4 5/8 x 3 1/2 x 1 1/8 in (11.7 x 8.9 x 2.8 cm)
B1200-13BK	4 5/8 x 3 1/2 x 1 1/8 in (11.7 x 8.9 x 2.8 cm)
B1200-15	4 1/2 x 4 1/2 x 1 1/4 in (11.4 x 11.4 x 2.8 cm)
B1200-15BK	4 1/2 x 1 1/2 x 1 1/4 in (11.4 x 11.4 x 2.8 cm)
B1200-17	6 x 4 x 1 1/4 in (15.2 x 10.2 x 2.8 cm)
B1200-19	6 x 4 x 2 in (15.2 x 10.2 x 5.1 cm)
B1206	6 compartment box (11.7 x 8.9 x 2.8 cm)