

HALO™

ADJUSTABLE PIPETTOR



The Halo™ Adjustable Pipettor

Halo™ Pipettors take a modern approach to classic liquid handling. As one of the lightest pipettors ever introduced, Halo Pipettors dramatically reduce the risk of repetitive stress injury and general shoulder, arm, and hand fatigue associated with daily pipetting. Reduced actuation forces and an ergonomics-first design make Halo Pipettors feel balanced and natural in hands of any size.

1

Expanded Volume Range

In addition to the traditional full spectrum volume range, the Halo pipettor is also available in an exclusive 1250uL model. Unlock the full potential of your extended length tips with this innovative new size.

2

DriftGuard™ Technology

The click-based volume adjustment wheel is designed to prevent volume drift. Even after hundreds of aspirations, the selected volume remains uncompromised.

3

Ultra Light-Weight

Halo Pipettors are one of the lightest adjustable pipettors ever made. At 40% lighter than the worldwide leader, Halo pipettors reduce the arm, wrist, and shoulder fatigue often associated with repetitive pipetting.

4


Calibration Simplified

Save hundreds on calibration services - Each Halo pipettor is equipped with a user enabled calibration port and adjustment key, which allow the user to perform in-lab calibrations quickly and accurately, using nothing more than an analytical balance.

3 year warranty

The Halo Advantage

Color Coded Tops
For Easy Identification



Key Features



The Lightest Pipettor Available
reduces repetitive stress and fatigue



4 Digit Micrometer
with DriftGuard™ protection



Exclusive 1250µL Size
unlock 25% more capacity from XL tips



Universal Tip Fit
compatible with all major brands



Pipettor Specifications

<u>Cat. No.</u>	<u>Volume Range</u>	<u>Accuracy</u>	<u>Reproducibility</u>	<u>Cat. No.</u>	<u>Volume Range</u>	<u>Accuracy</u>	<u>Reproducibility</u>
H6800-2	0.2-2µL	±12.0-2.5%	≤6.00-2.00%	H6800-200	20-200µL	±3.0-0.6%	≤1.00-0.15%
H6800-10	1-10µL	±2.5-1.0%	≤1.50-0.80%	H6800-1M	100-1000µL	±3.0-0.6%	≤0.60-0.20%
H6800-20	2-20µL	±3.0-0.9%	≤2.00-0.40%	H6800-1MXL	200-1250µL	±3.0-0.6%	≤0.60-0.20%
H6800-50	5-50µL	±2.0-0.6%	≤2.00-0.30%	H6800-5M	500-5000µL	±2.4-0.6%	≤0.60-0.15%
H6800-100	10-100µL	±3.0-0.8%	≤1.50-0.15%	H6800-10M	1000-10000µL	±3.0-0.6%	≤0.60-0.20%

Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Quantity</u>	<u>Price</u>
H6800-2	Halo™ Single Channel Pipettor, lightweight, 0.2-2µL	1/ea	
H6800-10	Halo™ Single Channel Pipettor, lightweight, 1-10µL	1/ea	
H6800-20	Halo™ Single Channel Pipettor, lightweight, 2-20µL	1/ea	
H6800-50	Halo™ Single Channel Pipettor, lightweight, 5-50µL	1/ea	
H6800-100	Halo™ Single Channel Pipettor, lightweight, 10-100µL	1/ea	
H6800-200	Halo™ Single Channel Pipettor, lightweight, 20-200µL	1/ea	
H6800-1M	Halo™ Single Channel Pipettor, lightweight, 100-1000µL	1/ea	
H6800-1MXL	Halo™ Single Channel Pipettor, lightweight, 200-1250µL	1/ea	
H6800-5M	Halo™ Single Channel Pipettor, lightweight, 0.5-5mL	1/ea	
H6800-10M	Halo™ Single Channel Pipettor, lightweight, 1-10mL	1/ea	
H6800-CR	Carousel Stand for 6 Halo pipettes	1/ea	
H6505	5mL Pipette Tips for Halo™ Pipettors 100/pack	100/pk	
H6510	10mL Pipette Tips for Halo™ Pipettors 100/pack	100/pk	
H6800-SK	Halo™ Starter kit includes 1 H6800-10, H6800-20, H6800-200, H6800-1MXL, and a carousel stand	1/kit	



MTC Bio - Proud to be a trusted supplier to life science



www.mtcbio.com



+1 (908) 279-0280



info@mtcbio.com