

# CleanPrep™ Dead Air Box - 110V

# MYSTAIRE®



Made in the USA



## GENERAL SPECIFICATIONS

<b>Construction Material</b>	<b>Work Surface</b>	Polypropylene
	<b>Shell</b>	Polypropylene
	<b>Front Sash</b>	Polycarbonate
	<b>Pipette Shelf</b>	Polycarbonate
<b>Controls</b>		ON/OFF Switch
<b>Lighting</b>		Fluorescent
<b>UV Light</b>		254-nm Shortwave UV Light
<b>UV Light Timer</b>		0-15 Minute UV Light Timer

## DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-DB24	24.0"	24.0"	29.5"	23.5"	20.0"	29.0"	22.0"	10.5"	23"
MY-DB36	36.0"	24.0"	29.5"	35.5"	20.0"	29.0"	34.0"	10.5"	23"
MY-DB48	48.0"	24.0"	29.5"	47.5"	20.0"	29.0"	46.0"	10.5"	23"

## TECHNICAL SPECIFICATIONS

Product #	Electrical		Unit Weight*
	Voltage	Power	
MY-DB24	110V	50 W	60 lbs.
MY-DB36	110V	50 W	75 lbs.
MY-DB48	110V	100 W	90 lbs.

Unit weight is an approximation\*

[www.mystaire.com](http://www.mystaire.com)

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: [contact@mystaire.com](mailto:contact@mystaire.com)