

Secure Dry™ Evidence Drying Cabinets - 110V

MYSTAIRE®



Made in the USA



GENERAL SPECIFICATIONS

Construction Material	Shell	Polypropylene
	Door	Polycarbonate Viewing Window Epoxy Powder-Coated Steel
	Filtration Housing	Polypropylene
	Shelves	Two - Stainless Steel
	Drying Rod	Stainless Steel
Controls		ON/OFF Switch Automatic Two-Speed Blower
Blowers		Vapor-Proof Fluorescent
Lighting		254-nm Shortwave UV Light 0-15 Minute Timer
Optional UV Light		Audible / Visual Alarm
Door Prefilter Monitoring		Audible / Visual Alarm
Filter Change Alarm		Audible / Visible Alarm

DIMENSIONS

Product #	External Dimensions			Internal Dimensions		
	Width	Depth	Unit Height	Width	Depth	Height
MY-FDR32	32"	23"	78"	31.25"	22.5"	55"
MY-FDR48	48"	23"	78"	47.25"	22.5"	55"
MY-FDR64	64"	23"	78"	63.5"	22.5"	55"
MY-FDR64D*	64"	23"	78"	31.4"	22.5"	55"
MY-FDR72	72"	23"	78"	71.5"	22.5"	55"
MY-FDR72D*	72"	23"	78"	35.4"	22.5"	55"

Internal unit width per chamber, unit has two chambers.*

TECHNICAL SPECIFICATIONS

Product #	Electrical		Volume of Filtered Air	Prefilters Required		HEPA Filters Required	Carbon Filters Required	Unit Weight**
	Voltage	Power		Door	Trace Evidence			
MY-FDR32	110V	240 W	250 CFM	1	1	1	1	210 lbs.
MY-FDR48	110V	240 W	250 CFM	1	1	1	1	300 lbs.
MY-FDR64	110V	260 W	250 CFM	2	2	2	2	390 lbs.
MY-FDR64D	110V	250 W	250 CFM	2	2	2	2	395 lbs.
MY-FDR72	110V	260 W	250 CFM	2	2	2	2	465 lbs.
MY-FDR72D	110V	250 W	250 CFM	2	2	2	2	475 lbs.

Unit weight is an approximation**

www.mystaire.com

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: contact@mystaire.com