


# Basic Ductless Workstation - 110V

# MYSTAIRE®



Made in the USA



| GENERAL SPECIFICATIONS |                    |  |
|------------------------|--------------------|--|
| Construction Material  | Base/Work Surface  | Optional, Removable Polypropylene Tray |
|                        | Shell              | Polycarbonate                          |
|                        | Front Sash         | Polycarbonate                          |
|                        | Filtration Housing | ABS                                    |
| Controls               |                    | User Controlled, Variable Air Volume   |
| Face Velocity          |                    | 60-100 LFM                             |
| Blowers                |                    | Variable Speed, Brushless-Sparkless    |

| DIMENSIONS |                     |       |        |                     |       |        |                         |         |           |
|------------|---------------------|-------|--------|---------------------|-------|--------|-------------------------|---------|-----------|
| Product #  | External Dimensions |       |        | Internal Dimensions |       |        | Sash Opening Dimensions |         |           |
|            | Width               | Depth | Height | Width               | Depth | Height | Width                   | Opening | Sash Fold |
| MY-2400    | 24"                 | 22.5" | 29.5"  | 23"                 | 22"   | 19"    | 23.5"                   | 8"      | 20.1"     |
| MY-3200    | 32"                 | 22.5" | 29.5"  | 31"                 | 22"   | 19"    | 31.5"                   | 8"      | 20.1"     |
| MY-4800    | 48"                 | 22.5" | 29.5"  | 47"                 | 22"   | 19"    | 47.5"                   | 8"      | 20.1"     |

| TECHNICAL SPECIFICATIONS |            |       |                        |                     |                       |              |
|--------------------------|------------|-------|------------------------|---------------------|-----------------------|--------------|
| Product #                | Electrical |       | Volume of Filtered Air | Prefilters Required | Main Filters Required | Unit Weight* |
|                          | Voltage    | Power |                        |                     |                       |              |
| MY-2400                  | 110 V      | 110 W | 175 CFM                | 1                   | 1                     | 45 lbs.      |
| MY-3200                  | 110 V      | 220 W | 480 CFM                | 1                   | 1                     | 80 lbs.      |
| MY-4800                  | 110 V      | 200 W | 350 CFM                | 2                   | 2                     | 120 lbs.     |

\*Unit weight is an approximation

[www.mystaire.com](http://www.mystaire.com)

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: [contact@mystaire.com](mailto:contact@mystaire.com)