


VEE GEE®

Crucibles

VEE GEE®, made by  Porcelain High Form Crucibles are manufactured with the highest standards. Crucibles are glazed inside and out except for outside bottom surface and rim. Crucible Covers with ring handle are glazed inside and out except for rim.



Specifications

Crucibles w/out Covers

Cat. No.	Capacity (mL)	Dimensions Top O.D. x Height	Pack Qty
52120-0301	10	30 x 25 mm	5
52120-0351	12	35 x 28 mm	5
52120-0451	30	45 x 36 mm	5
52120-0501	45	50 x 40 mm	5
52120-0601	80	60 x 48mm	5
52120-0701	120	70 x 56 mm	5



Covers Only

Cat. No.	Fits Diameter	Pack Qty
52100-0301	30 mm	5
52100-0351	35 mm	5
52100-0451	45 mm	5
52100-0501	50 mm	5
52100-0601	60 mm	5
52100-0701	70 mm	5



VEE GEE

VEE GEE Scientific, Inc.

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084