

# Tygon® A-60-F Tubing

## Food and Beverage

### A-60-F Tubing

Tubing will not crack or deteriorate, even in extreme temperatures. Featuring high flexural fatigue resistance, tubing is perfect for use in peristaltic pumps often found in food and beverage dispensing equipment. It resists kinks, retains its shape, and remains flexible for quick and easy installation.

This high-performance tubing withstands harsh sanitizing solutions and exhibits exceptional acid and alkali resistance.

Easily clean tubing using clean-in-place or steam-in-place cleaning and sterilization systems or autoclave up to five cycle times without affecting the service life.

### APPLICATIONS

- General-purpose, laboratory, and single-use applications

### FEATURES/BENEFITS

- Free technical application support available to assist with tubing selection
- Excellent flexural fatigue resistance—great for use in peristaltic pumps
- Steam clean in place to eliminate frequent tubing replacement
- Chemically compatible with a variety of cleaning solutions



### CERTIFICATION

Cole-Parmer will provide a certificate of compliance FREE at your request where applicable.

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	86.21.5109.9909
France	+33 (0) 1486 37800

Germany	+49 (0) 9377 92030
India	+1.800.266.1244
Italy	+39 (0) 1313 89513
UK	+44 (0) 1480.272279
All others	+1.847.549.7600

## Specification and Ordering Table

Dimensions: in. (mm)			Max pressure: PSI (bar)	Item number	Roll length (ft)
Tubing ID	Tubing OD	Wall			
5/8"	7/8"	0.0625 (1.6)	34 (2.3)	50106-08	50
3/8"	1/2"	0.0625 (1.6)	19 (1.3)	50106-05	50
1/2"	3/4"	0.0625 (1.6)	13 (0.9)	50106-07	50
1/4"	1/2"	0.0625 (1.6)	10 (0.7)	50106-03	50
1/16"	3/16"	0.125 (3.2)	19 (1.37)	50105-99	50
3/4"	1"	0.0625 (1.6)	8 (0.55)	50106-09	50
1/8"	1/4"	0.0625 (1.6)	7 (0.48)	50106-00	50
5/16"	7/16"	0.125 (3.2)	13 (0.9)	50106-04	50
3/8"	5/8"	0.125 (3.2)	10 (0.7)	50106-06	50
1/4"	3/8"	0.125 (3.2)	8 (0.55)	50106-02	50
3/16"	5/16"	0.125 (3.2)	7 (0.48)	50106-01	50



[sales@novatech-usa.com](mailto:sales@novatech-usa.com)

[www.novatech-usa.com](http://www.novatech-usa.com)

Tel: (866) 433-6682

Fax: (866) 433-6684

Tel: (281) 359-8538

Fax: (281) 359-0084

**Cole-Parmer®**  
an antylia scientific company