

United Scientific™ Float and Spin Magnet Demonstrator



This attention-getting demonstration of magnetic repulsion uses a clear acrylic plate that keeps a spinning axle oriented above magnets built into the base of the apparatus. With a simple twist of your fingers, you can start this levitating magnet spinning and kick start a classroom discussion. The black and white coloring makes the demonstration easy to see and helps track the speed of the rotation.

The axle measures 5" long Dimensions: 2-1/2" x 5" x 3"

- **Easy to see black and white pattern**
- **Smooth rotation**
- **No batteries required**



Black and white pattern makes the rotation easy to see.

Item #	Description
FLSPMAG	United Scientific™ Float and Spin Magnet Demonstrator

⚠ WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov.



**United
Scientific®**

1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556