

United Scientific™ Glass Erlenmeyer Flask Set of 3



An ideal starter set for any laboratory, this set contains three narrow-mouth, durable borosilicate Erlenmeyer flasks in sizes of 50 mL, 150 mL, and 250 mL. Borosilicate glass flasks are designed to give the best combination of thermal shock resistance and physical strength. Features heavy-duty rims and graduations in durable white enamel.

- **Excellent thermal shock resistance**
- **Durable white enamel markings**
- **Heavy-duty rims**
- **Autoclavable**
- **Tolerances to $\pm 5\%$**

Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Tolerance (\pm mL)	Stopper Size
50	20 - 50	10	$\pm 5\%$	3
150	75 - 150	25	$\pm 5\%$	5
250	50 - 200	25	$\pm 5\%$	6.5

Item #	Description
FGSET3	United Scientific™ Glass Erlenmeyer Flask Set of 3



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556