

United Scientific™ Laser Pointers, Class 2



Economical laser pointers are useful for many optical experiments. Ideal for middle school and high school classrooms. Lasers operate on a push-button switch, are housed in a rubberized aluminum pen-style case, and include a pocket clip.

The 6" long x 0.5" diameter Class 2 Red laser pointer has a maximum output of 1mW and a wavelength of 650nm (+/- 10nm). Range of beam is 75 meters. Operates on two "AAA" batteries (included).

The 6" long x 0.5" diameter Class 2 Green laser pointer has a maximum output of 1mW and a wavelength of 532nm (+/- 10nm). Range of beam is 125 meters. Operates on two "AAA" batteries (included).



RLPCLS2 Red Laser Pointer, Class 2

- **Maximum output of less than 1 mW**
- **Push-button switch**
- **Red and green versions available**



GLPCLS2 Green Laser Pointer, Class 2

Item #	Description
RLPCLS2	United Scientific™ Red Laser Pointer, Class 2
GLPCLS2	United Scientific™ Green Laser Pointer, Class 2

 **WARNING:** Cancer and Reproductive Harm - www.P65Warnings.ca.gov.



**United
Scientific®**

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com
Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084