



Variable Volume Fully Autoclavable Micropipettes

Our Variable Volume Micropipettes feature built-in tip ejectors and autoclavable tip cones. They are offered in nine volume ranges, as detailed below. They work on a click stop digital system, are easy to calibrate and maintain, and easy to disassemble for autoclaving. Manufactured in an ISO 9001:2008 certified facility, these pipettes are highly accurate. Each pipette is individually calibrated according to ISO 8655 standards. A calibration certificate is provided with each pipette.



| Item No. | Volume Range (µl) | Increment (µl) | Volume (µl) | Accuracy (± %) | Precision (± %) | Compatible Tips |
|-----------|-------------------|----------------|-------------|----------------|-----------------|-----------------|
| PRV-2 | 0.2-2.0 µl | 0.01 | 0.2 | 12.0 | 6.0 | 20 µl |
| | | | 1.0 | 2.5 | 1.5 | |
| | | | 2.0 | 1.5 | 0.7 | |
| PRV-10 | 0.5-10 µl | 0.1 | 0.5 | 5.0 | 2.8 | 200 µl |
| | | | 5.0 | 1.5 | 0.8 | |
| | | | 10.0 | 1.0 | 0.4 | |
| PRV-20 | 2-20 µl | 0.1 | 2.0 | 5.0 | 1.5 | 200 µl |
| | | | 10.0 | 1.0 | 0.6 | |
| | | | 20.0 | 1.0 | 0.3 | |
| PRV-50 | 5-50 µl | 1.0 | 5.0 | 2.5 | 1.5 | 200 µl |
| | | | 25.0 | 0.7 | 0.3 | |
| | | | 50.0 | 0.7 | 0.3 | |
| PRV-100 | 10-100 µl | 1.0 | 10.0 | 2.5 | 0.7 | 200 µl |
| | | | 50.0 | 0.8 | 0.3 | |
| | | | 100.0 | 0.8 | 0.15 | |
| PRV-200 | 20-200 µl | 1.0 | 20.0 | 2.5 | 0.7 | 200 µl |
| | | | 100.0 | 1.0 | 0.3 | |
| | | | 200.0 | 0.6 | 0.2 | |
| PRV-1000 | 100-1000 µl | 10.0 | 100.0 | 1.6 | 0.3 | 1000 µl |
| | | | 500.0 | 0.6 | 0.2 | |
| | | | 1000.0 | 0.5 | 0.2 | |
| PRV-5000 | 0.5-5 ml | 100.0 | 500.0 | 1.5 | 0.3 | 5 ml |
| | | | 2500.0 | 0.6 | 0.2 | |
| | | | 5000.0 | 0.6 | 0.2 | |
| PRV-10000 | 2-10 ml | 100.0 | 1000.0 | 1.5 | 0.3 | 10 ml |
| | | | 5000.0 | 0.6 | 0.2 | |
| | | | 10000.0 | 0.6 | 0.2 | |



Easy Volume Adjustment:
Easy volume adjustment by simply turning the plunger. Plunger does not snag gloves.



Volume Setting:
Soft click sound at every increment ensures perfect volume setting and no accidental volume change. It also facilitates single handed operation.



Digital Display:
Four-digit display with subdivisions provides small volume increments.



Color Coding:
Different volume range micropipette dials are color-coded for easy identification.



Large Grip:
Specially designed large grip is provided for comfort during frequent use.



Tip Ejector:
Built-in streamlined tip ejector facilitates easy tip ejection and access to narrow necked bottles and tubes.



User Adjustment:
Re-calibration can be performed.

Contact us for information regarding private labeling.