

United Scientific™ PinPoint™ Pipette Controller



Ultra light and powerful pipette controller is designed for comfort as well as accuracy. Intuitive 6 level speed setting and sensitive triggers provide a wide range of aspirating/dispense speeds. Pipette controller is ergonomically designed for extended usage; the lithium-ion battery charges quickly and can provide up to 8 hours use.

- **Compatible with glass or plastic serological pipettes (.1 mL - 100 mL)**
- **LCD indicates battery charge and 6 speed levels**
- **Fills a 25 mL pipette in less than 5 seconds**
- **Lithium-ion battery charges fully within 3 hours and can provide up to 8 hours use***
- **Accessories include: desk stand, wall mount, extra filter, and AC adapter with USB connection**
- **Meets the requirements of the applicable EU directives, UK regulations, UL, and CAN / CSA**



United Scientific™ PinPoint™ Pipette Controller

LCD indicates battery charge level and aspirate / dispense speed settings

Item #	Description	Quantity
UNPINPNT01	United Scientific™ PinPoint™ Pipette Controller	1
UNPINPNT02	United Scientific™ PinPoint™ Silicone Adapter	1
UNPINPNT03	United Scientific™ PinPoint™ Nosecone	1
UNPINPNT04	United Scientific™ PinPoint™ .45 µm Filters	5
UNPINPNT05	United Scientific™ PinPoint™ .22 µm Filters	5

* Charge time shown is after fully charging for first time use. Battery usage can vary based on continuous or intermittent use.



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682 Fax: (866) 433-6684
Tel: (281) 359-8538 Fax: (281) 359-0084

United Scientific™ PinPoint™ Pipette Controller



What's in the box


PinPoint Pipette Controller with nosecone / silicone adapter / .45 µm filter, desk stand, wall mount, extra .45 µm filter, AC adapter with USB connection, and a user manual.



Compatible with glass or plastic
serological pipettes (.1 mL - 100 mL)



Desk stand



Wall mount

Parts also available separately:



.45 µm and
.22 µm filters



Nosecone



Silicone adapter



PinPoint can be oriented in two different
positions on the desk stand.



PinPoint can also be easily
holstered into the wall mount.

Recommended accessory: Sterile, Disposable Serological Pipettes

- Individually wrapped
- Gamma sterilized
- Barrier plug to prevent overflow
- Sizes available from 1 mL to 50 mL

See our product offering [here](#)



1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556

novatech
INTERNATIONAL

sales@novatech-usa.com
www.novatech-usa.com

Tel: (866) 433-6682
Tel: (281) 359-8538

Fax: (866) 433-6684
Fax: (281) 359-0084