

# United Scientific™ SteriTrue™ PS Sampling Scoops

## Specifications

### Description:

- SteriTrue™ Sterile Sampling Scoop, Polystyrene, 125 mL (4 oz)

### Features:


- Ergonomically designed with voice-of-customer input
- Flat bottom surface for stability during weighing
- Gamma sterilized, individually packaged
- Lot stamped for traceability
- PS material meets FDA guidelines

### Physical Properties:

- Material: High-impact polystyrene
- Meets FDA criteria in 21 CFR 177.1520 (2017)
- Heavy Metals free, BPA & Phthalate free, BSE/TSE compliant
- Packaging: sealed polyethylene bag
- Sterilization: Gamma radiation
- Expiration date: 2 yrs from sterilization date
- Overall Length: 221.30mm (7.71 inches)
- Handle Length: 126.50mm (4.98 inches)
- Maximum Width: 54.25mm (2.14 inches)
- Weight: 1.05oz (29.8g)
- Quantity per pack: 100/Pk and 10/Pk
- Country of origin: India



Item #	Description	Quantity per pack
UNST101	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 125 mL (4 oz)	100
UNST201	SteriTrue™ Sterile Sampling Scoops, Polystyrene, 125 mL (4 oz)	10

 **WARNING:** This product can expose you to chemicals including Styrene, which is [are] known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).